

# New! Stander Size 1



**NEW**



# Stander

**It's finally here** – the Stander that has it all. We're thrilled to introduce the first size of a total redesign of our Supine and Prone Stander families. This new Stander can be configured as supine, prone or multi-position.

All models allow hip abduction up to 30° to promote hip joint development. Independent adjustments for leg length can accommodate length discrepancies. A broad range of positioning components allow users to remain comfortably upright and active for longer. The frame can be raised and lowered to position a child at eye level with peers, or in easy reach of activities.

All this function is contained within an elegant and approachable design that centers the child while opening worlds of possibility.



Find letters of medical necessity and studies on the therapeutic benefits of standing at:  
[www.rifton.com/newstander](http://www.rifton.com/newstander)

*"The* ability to fully customize joint positioning through the lower extremities for maximal comfort and alignment is a game changer! I have a student who only tolerated five to ten minutes of static standing in a less adjustable stander. When we tried the new Rifton Stander the student was happy and comfortable for nearly an hour, and we only got him out because we ran out of time."*"*

*Linda K. Miller, PT, DPT  
MOVE International Trainer  
Prospect Center, NY*

●..... **Supine Stander configuration** (S410)

# Stander features



**Supine Stander**



The Stander has abduction up to 30°.



The Stander has 0°–90° gas-spring-assisted tilt adjustment.



### Prone Stander

# Multi-Position Stander



The Stander can be configured in prone or supine. If the multi-position configuration is ordered, the conversion from one to the other is quick and easy.

## Quick and easy conversion



**Prone**

- 1 **Add/Remove** the upper trunk support.
- 2 **Rotate** the tray.
- 3 **Swap** the supine hip strap and prone pelvic harness.
- 4 **Add/Remove** the supine knee straps.
- 5 **Turn** the sandals 180°.



**Supine**



## Hip abduction

The Stander supports abduction up to 30°. Research has shown that weight bearing in an abducted position can improve hip joint health and development during childhood, particularly for children at risk of hip displacement.



When the frame is lowered, users can interact at eye level with their peers.



With the frame raised, users can access activities and environments.





Sit-to-stand transfers are easy in the prone configuration.



Pads and supports swing away for easy lateral transfers in the supine configuration.

# Components



## Knee cuffs (prone only)

The knee cuffs provide lateral stabilization and rotational control and prevent knee hyperextension.



## Head laterals (supine only)

The head laterals mount on the headrest to provide lateral head support. They adjust independently in height and width.



## Butterfly harness (supine only)

The butterfly harness can be used instead of the trunk strap for greater upper trunk stabilization.



Standard



Rotation control

## Hip strap (supine only)

This strap supports the user at the level of the hips. The rotation control version provides rotational stabilization.



Standard



Rotation control

## Trunk strap

The trunk strap supports the user's chest. The rotation control version wraps around the user and provides rotational stabilization.





### Tray

The tray can be used in prone and supine configuration. Its position and angle adjust easily and it swings away for transfer.




### Additional trunk laterals

The laterals attach at chest or hip level and swing away for transfer. The pads rotate for individual height adjustment.




### Medial/lateral thigh supports

The medial/lateral thigh supports provide abduction or adduction. The mounting bracket is off-center, so the pad has two possible height positions depending which side it is installed on.



# Stander dimensions

User dimensions (inches)	S410 Size 1
Height	28–42
 <p><b>Key user dimension: height</b> Select the appropriate stander by the user's overall height. Choose the model that allows for growth.</p> <p><b>Important:</b> User's weight must not exceed the maximum working load.</p>	
Item dimensions (inches)	S410 Size 1
Base length x width	32½ x 21½
Height of horizontal body support above floor	22–28
Sandal to knee height	6¾–11
Knee to hip height	6¼–10½
Sandal to top of headrest	31½–43
Width between trunk laterals	5½–10
Width between medial/lateral thigh supports	2½–4
Basic item weight PRONE (lb)	47
Basic item weight SUPINE (lb)	52
Max. working load (lb)	75



**Range of adjustment**

# 2024 order form

Use dimension chart to select appropriate size.

			S410 size 1		
<b>* Frame</b>			S411	<input type="checkbox"/>	
<b>* Color</b>			<b>CHOOSE</b>	Red <input type="checkbox"/>	
				Blue <input type="checkbox"/>	
<b>* CHOOSE CONFIGURATION</b>	<b>Prone</b> (with laterals, pelvic harness)		S412	<input type="checkbox"/>	
	HCPCS code: E0638	<b>Prone knee cuffs</b> (pair)	S485	<input type="checkbox"/>	
	<b>Supine</b> (with laterals, knee strap)		S413	<input type="checkbox"/>	
	HCPCS code: E0638	<b>Head laterals</b> (pair)	S497	<input type="checkbox"/>	
	<b>Butterfly harness</b>		S456	<input type="checkbox"/>	
	<b>* Hip strap</b>	<b>CHOOSE</b>	Standard	S472	<input type="checkbox"/>
			Rotation control	S476	<input type="checkbox"/>
	<b>Multi-Position</b> (with laterals, knee strap, pelvic harness)		S414	<input type="checkbox"/>	
	HCPCS code: E0641	<b>Head laterals</b>	S497	<input type="checkbox"/>	
	<b>Butterfly harness</b>		S456	<input type="checkbox"/>	
<b>* Hip strap</b>	<b>CHOOSE</b>	Standard	S472	<input type="checkbox"/>	
		Rotation control	S476	<input type="checkbox"/>	
<b>Prone knee cuffs</b> (pair)		S485	<input type="checkbox"/>		
<b>* Trunk strap</b>			<b>CHOOSE</b>	Standard <input type="checkbox"/>	
				Rotation control <input type="checkbox"/>	
<b>* Leg type</b>			<b>CHOOSE</b>	Standard <input type="checkbox"/>	
				Advanced <input type="checkbox"/>	
<b>Tray</b>			S491	<input type="checkbox"/>	
<b>Additional trunk laterals</b> (pair)			S494	<input type="checkbox"/>	
<b>Medial/lateral thigh supports</b> (pair)			S482	<input type="checkbox"/>	
<b>Additional medial/lateral thigh supports</b> (pair)			S482	<input type="checkbox"/>	
<b>Total</b>			<b>S410</b>	<b>\$ _____</b>	

**\* Required**

The HCPCS codes indicated above are suggestions only, based on knowledge of our products and the HCPCS definitions established by CMS. It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. Actual coding and coverage varies by payer; some may accept or require alternative HCPCS codes, including miscellaneous codes, to ensure access for their beneficiaries.



Rifton Equipment  
PO Box 260  
Rifton NY 12471-0260



▶ Visit [www.rifton.com](http://www.rifton.com) to get a quote! 800.571.8198