



Convaid catalogue


Cruiser, EZ Rider & Rodeo

 **etac**[®]



Seating

Convaid



User height
600–1800 mm



Max user mass
114 kg



Environment
Indoor & outdoor

Cruiser & Cruiser Scout

The Cruiser features 30° of fixed tilt for improved postural control. Seat depth adjustment and self-tensioning seating for increased sitting tolerance, the durable and convenient Cruiser is the perfect solution for an adventurous family.

Highlights about Cruiser



Customized growth capabilities

Adjustable seat depth extends chair life for a growing child.



30° fixed tilt

Increased tilt angle improves postural control through upper body positioning.



Foldable

The Cruiser is folded quickly and fits into almost any car.



Self-Tensio®

Built in self-tensioning seating system increases sitting tolerance.



Standard, Cruiser

30° fixed tilt, removable height adjustable swing-away footplates, seat depth growth of 130 mm, cordura upholstery, Self-tension seat, 2-point positioning belt, closure Strap.

Standard, Cruiser Scout

30° fixed tilt, removable height adjustable swing-away footplates, seat depth growth of 130 mm, heavy-duty reinforced cordura, Self-tension seat, anatomic back, flip up arm rest, 2-point positioning belt, closure Strap.

Options, Cruiser

Sizes:

- CX10
- CX12
- CX14
- CX16
- CX18

Brakes

- Foot brake
- Foot brake and attendant hand brakes (not for CX 18)

Tyres

- Solid
- Air
- Solid w/transport fitting
- Air w/transport fitting

Options, Cruiser Scout

Sizes

- CX12
- CX14
- CX16
- CX18

Tyres

- Air
- Air, w/transport fitting

Cushion colours



Royal blue



Cherry red



Panther black*



Pink



Purple

Frame colour



Silver

*Standard in AU














Cruiser Scout 16"



Cruiser standard 18"

Measurements

		CX 10	CX 12	CX 14	CX 16	CX 18
	Effective seat width	255 mm	305 mm	355 mm	405 mm	455 mm
	Effective seat depth	152-299 mm	200-330 mm	254-380 mm	279-406 mm	356-533 mm
	Backrest height	530 mm	560 mm	640 mm	710 mm	640 mm
	Mass device	12 kg	12 kg	13 kg	14kg	14 kg
	Max user mass regular/transit	34/30 kg	34/30 kg	45,5/45,5 kg	77/77 kg	114/77 kg
	Overall width	560 mm	600 mm	650 mm	700 mm	724 mm
	Length	1220 mm	1120 mm	1245 mm	1320 mm	1400 mm
	Overall height	940 mm	970 mm	940 mm	1020 mm	970 mm
	Folded length	1190 mm	1220 mm	1190 mm	1350 mm	1270 mm
	Folded width	380 mm	380 mm	410 mm	410 mm	380 mm
	Folded height	430 mm	430 mm	430 mm	430 mm	430 mm

Further measurements at our website

Nice to know

The Cruiser is lightweight and has an easy three dimensional fold, which mean it requires minimal storage space and is easy to take along.

The individual swing away footplates allows for standing transfer.



Accessories



Folding top with windows

CX10: C904736-xx
CX12: C904737-xx
CX14: C904738-xx
CX16: C904739-xx
CX18: C904740-xx




Rain cover

CX10-12: C904939
CX14: C904940
CX16-18: C905283



Headrest extension

Curved cordura
CX10: C903332-xx
CX12: C902863-xx
CX14: C900484-xx
CX16: C900138-xx
CX18: C902619-xx



Headwings

Padded:
CX10: C903331
CX12: C902862
CX14: C900681
CX16: C900142
CX18: C902618



Occi headwings

C90527x



3-pt positioning belt

C904017
With depth adjustment croch strap



H-harness with padded covers

C903988
3-pt position belt required



Full support torso vest

CX10: C903330
CX12: C902861
CX14: C900680
CX16: C900339
CX18: C902616
3 pt. position belt required



Soft lateral support

CX10: C903328
CX12: C902859
CX14: C900678
CX16: C900140
CX18: C902614



Reducer seat insert

Seat width 5 cm smaller, seat depth 2,5 cm smaller
CX10-12: C904211
CX14-18: C904110



Upper extremity support surface

Clear tray
CX10: C903078
CX12: C903385
CX14: C903386
CX16: C903387
CX18: C903388



Medial thigh support-abductor

With position/align cushion
CX10: C903326
CX12: C902857
CX14: C900480
CX16: C900137
CX18: C902607

For more accessories - visit our website



Lateral thigh support-abductor

with position/align cushion
 CX10: C903327
 CX12: C902858
 CX14: C900482
 CX16: C900188
 CX18: C902608



Calf panel

C904580



Transport fittings

(Always together with Headrest Extension)
 CX10-16: C997125
 CX18: C997141
 Scout model:
 CX12-16: C997219
 CX18: C997220



Saddle bag

C904094




Foot positioners

Small (18-25 cm): C903415
 Medium (20-28 cm): C903416
 Large (23-30 cm): C903417



Angle adjustable footplates

CX10: C903445
 CX12-14: C903268
 CX16-18: C903269




Footplate securement strap

C903584



Under seat storage basket

CX10-12: C904046
 CX14-16: C904047
 CX18: C904048



Anti-tip

CX10-12: C904520
 CX14-16: C904540



Curb tipper

CX10-16: C904128
 CX18: C903916
 Scout: C903916



Utility bag

C904531



Wheelchair travel bag

C902022



Seating

Convaid



User height
1000–1800 mm



Max user mass
114 kg



Environment
Indoor & outdoor

EZ Rider

The Convaid EZ Rider is a pediatric buggy with 10° fixed tilt. It is lightweight and has compact folding. The EZ Rider has a huge range of accessories to meet the need of each individual child.

Highlights about EZ Rider




Customized growth capabilities

Adjustable seat depth extends chair life for a growing child.



10° fixed tilt

Provides upright posture to enhance functional positioning.



Upholstery & colours

Removable, washable Cordura available in various colours.



Foldable

The EZ Rider is folded quickly and fits into almost any car.



Standard

10° fixed tilt, removable height adjustable swing-away footplates, seat depth growth of 130 mm, cordura upholstery, Self-tension seat, closure bar, belt.

Options

Sizes:

- EZ12
- EZ14
- EZ16
- EZ18

Tyres

- Solid
- Air
- Solid w/transp. fitting
- Air w/transp. fitting

Cushion colours



Blue

*Standard in AU



Red



Black*



Pink



Purple

Frame colour



Silver




EZ Rider Standard, size 12" - 16"



EZ Rider Standard, size 18"

Measurements

		EZ 12	EZ 14	EZ 16	EZ 18
	Effective seat width	305 mm	355 mm	405 mm	455 mm
	Effective seat depth	230-360 mm	280-410 mm	305-430 mm	370-495 mm
	Backrest height	530 mm	610 mm	660 mm	680 mm
	Mass device	11 kg	14 kg	14kg	13 kg
	Max user mass regular/transit	34/30 kg	45,5/45,5 kg	77/77 kg	114/77 kg
	Overall width	546 mm	548 mm	660 mm	660 mm
	Length	915 mm	950 mm	1105 mm	1220 mm
	Overall height	965 mm	980 mm	1015 mm	950 mm
	Folded length	1130 mm	1170 mm	1220 mm	1160 mm
	Folded width	360 mm	340 mm	360 mm	370 mm
	Folded height	445 mm	460 mm	430 mm	390 mm

Further measurements at our website

Nice to know

EZ Rider is lightweight, portable and folds compactly. It is suitable for children in need of a functional upright seating position, where the focus is on allowing the child to be active and aware of the surroundings.



Accessories



Folding top with windows

EZ12: C904737-xx
EZ14: C904741-xx
EZ16: C904742-xx
EZ18: C904743-xx



Headrest cover

EZ12: C900853
EZ14: C901210
EZ16: C901205
EZ18: C901211



Rain cover

12: C904939
14: C904940
16/18: 905253



Headrest extension

Curved cordura
EZ12: C902863-xx
EZ14: C900484-xx
EZ16: C900138-xx
EZ18: C902619-xx



Headwings

Padded:
EZ12: C900849
EZ14: C900682
EZ16: C900683
EZ18: C900684



Occi headwings

C90527x



H-harness with padded covers

C903988
3-pt position belt required




3-pt positioning belt

C904017
With depth adjustment croth strap



Full support torso vest

EZ12: C900848
EZ14: C900338
EZ16: C901020
EZ18: C900387



Soft lateral support

Single
EZ12: C900846
EZ14: C900327
EZ16: C901018
EZ18: C900386



Reducer seat insert

Seat width 5 cm smaller
Seat depth 2,5 cm smaller
EZ12: C904211
EZ14-18: C904110



Height adjustable flip-up armrests

EZ12-14: C904523
EZ16-18: C904524

For more accessories - visit our website



Upper extremity support surface
Clear tray
EZ12: C904525
EZ14: C904526
EZ16: C904527
EZ18: C904528




Upper extremity support surface
Foam tray
EZ12: C903669
EZ14-16: C903670
EZ18: C903671



Medial thigh support-abductor
with position/align cushion
EZ12: C900843
EZ14: C900329
EZ16: C901024
EZ18: C900388



Lateral thigh support-abductor
with position/align cushion
EZ12: C900844
EZ14: C900330
EZ16: C901026
EZ18: C900389



Calf panel
C904580




Transport fittings
(Always together with Headrest Extension)
C997267




Foot positioners
Small (18-25 cm): C903415
Medium (20-28 cm): C903416
Large (23-30 cm): C903417



Angle adjustable footplates
EZ12/14: C903268
EZ16/18: C903269



Footplate securement strap
C903584



Anti-tip
EZ12-14: C900377
EZ16-18: C900378



Utility bag
C904531



Wheelchair travel bag
C902022



Seating

Convaid



User height
700–1700 mm



Max user mass
77 kg



Environment
Indoor & outdoor

Rodeo

The Rodeo is a buggy with 45° tilt-in-space and back recline. This buggy provides growth options through its adjustable seat depth feature. Rodeo has a huge range of accessories to meet the need of each individual child.

Highlights about Rodeo



Adjustable tilt and back recline

Provides pressure relief and support.



Customisable growth capabilities

Adjustable seat depth extends chair life for a growing child.



Lightweight and durable

Lightweight frame combined with high Max user mass - up to 77 kg.



Individual positioning

Accommodates functional or rest positions and increased sitting tolerance for maximum comfort.



Standard

5° to 45° adjustable tilt, 90° to 110° back recline, elevating legrest - one-piece, angle adjustable footplate, closure Strap.

Options

Sizes:

- RD12
- RD14
- RD16

Brake

- Foot brake
- Attendant hand brakes

Tyres

- Solid
- Air
- Solid w/transp. fitting
- Air w/transp. fitting

Transport Fittings

- With
- Without

Cushion colours



Blue



Red



Black*



Pink



Purple

*Standard in AU

Frame colour



Silver




Standard Rodeo RD14



Folded Rodeo

Measurements

		RD 12	RD 14	RD 16
	Effective seat width	305 mm	355 mm	405 mm
	Effective seat depth	230-360 mm	280-410 mm	330-460 mm
	Backrest height	560 mm	610 mm	635 mm
	Mass device	15 kg	16 kg	17 kg
	Max user mass regular/transit	34/30 kg	46/46 kg	77/64 kg
	Overall width	580 mm	630 mm	680 mm
	Length	890 mm	890 mm	910 mm
	Overall height	1090 mm	1120 mm	1140 mm
	Folded length	860 mm	860 mm	1015 mm
	Folded width	585 mm	635 mm	685 mm
	Folded height	430 mm	460 mm	460 mm

Further measurements at our website

Nice to know

Rodeo is suitable for walks in uneven terrain and nature.

It is easy to adjust and mount different accessories which makes it possible to use for children with different needs.



Accessories



Folding top with windows

12: C904745-xx
14: C904746-xx
16: C904747-xx



Headrest cover

C902679-xx
C902659-xx
C903152-xx



Rain cover

12: C904939
14: C904940
16: C905253



Headrest Extension

12: C902863-xx
14: C900484-xx
16: C900138-xx



Headwings

Padded:
12: C902673
14: C902653
16: C903147



Occi headwings

C90527x



H-harness with padded covers

With padded covers
C903988



3-pt positioning belt

With adjustable crotch strap
C904017



Full torso vest

12: C903261
14: C903260
16: C903259




Reducer seat insert

Seat width 5 cm smaller, seat depth 2,5 cm smaller.
12: C904211
14: C904110



Swing away laterals - contoured

X-small: C902980
Small: C902981
Medium: C902982
Large: C902983



Solid/Mesh side panels

12: C902672/ C902680
14: C902652/ C902681
16: C903145/ C903146

For more accessories - visit our website



**Height adjustable
flip-up armrest**

C904346



**Upper extremity support surface
(clear tray)**

Flip-up arm rest required

12: C903405

14: C903401

16: C903402



Upper extremity support tray

Foam tray

12: C903669

14/16: C903670



Medial thigh support-abductor

With position/align cushion

C902667

C902647

C903140




Lateral thigh support-abductor

With position/align cushion

C902668

C902648


C903141



Transport fittings

Always together with Headrest Extension

C997134



Legrest extension

12: C904095

14: C904096

16: C904097



Rear anti-tip tubes (pair)

C904649



Under seat storage basket

12: C904529


14/16: C904530



Medical necessity storage

12: C904476

14/16: C904477



Oxygen tank bag

C902870



Utility bag

C904531



© 2021 | Brochure, Gohard - AU | 2021 | Reserving our rights for printer's error or discontinued products.
Please find an updated edition on our website © R32 A/S All rights reserved.

Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our heart lies in the solutions that optimise quality of life for the individual, their family and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com

etac

Etac
PO Box 284
Botany NSW 1455
Australia
Tel +61 (02) 8213 6666
Fax +61 (02) 8213 6655
www.etac-au.com
au.sales@etac.com

etac[®]
Creating Possibilities