

In-situ Deluxe Leg Sling HS

S16

- Designed to provide full body support on transfer
- Adjustable integral head support with removable stiffeners
- Shaped leg sections provide full femur & hip support
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Detachable colour coded loops
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)



*nb: Sling shown in image is an In-Situ Deluxe Leg Sling HS with fixed Silva-Safety Slots**

Product Codes & Specifications:



	Silva-Superfine PLUS®	Silva-Superfine® with COOLMAX®	Silva-Superfine® with COOLMAX® Double Layer	A	B	C	D
XXXS	S16FXXXSM	S16SFCXXXSM	S16SCDLXXXS	56 - 66cm	50 - 55cm	56cm	60cm
XXS	S16FXSM	S16SFCXXSM	S16SCDLXXS	66 - 76cm	56 - 62cm	64cm	65cm
XS	S16FXSM	S16SFCXSM	S16CDLXS	76 - 86cm	60 - 66cm	73cm	70cm
S	S16FSM	S16SFCSM	S16SCDLS	86 - 96cm	66 - 74cm	80cm	78cm
M	S16FMED	S16SFCMED	S16SCDLM	96 - 108cm	74 - 80cm	90cm	86cm
L	S16FLRG	-	S16SCDLL	108 - 130cm	80 - 88cm	97cm	94cm
XL	S16FXLRG	-	S16SCDLXL	130 - 144cm	88 - 100cm	114cm	103cm
XXL	S16FXLRG	-	S16SCDLXXL	144 - 160cm	96 - 104cm	120cm	107cm

