

In-situ Straight Leg Sling HS

S15

- Designed to provide full body support on transfer
- Adjustable integral head support with removable stiffeners
- Easy to apply with low-friction fabric to leg areas and head support for ease of application
- Soft binding as standard
- Detachable colour coded loops
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 71°C/159°F
- Tumble dry low heat (max 50°C)



nb: Sling shown in image is an In-Situ Straight Leg HS sling with fixed Silva-Safety Slots*

Product Codes & Specifications:



	Silva-Superfine PLUS®	Silva-Superfine® with COOLMAX®	Silva-Superfine® with COOLMAX® Double Layer	A	B	C	D
XXXS	S15FXXXSM	S15SFCXXXSM	S15SFCDLXXXS	56 - 66cm	50 - 55cm	56cm	60cm
XXS	S15FXXSM	S15SFCXXSM	S15SFCDLXXS	66 - 76cm	56 - 62cm	64cm	65cm
XS	S15FXSM	S15SFCXSM	S15SFCDLXSM	76 - 86cm	60 - 66cm	74cm	70cm
S	S15FSM	S15SFCSM	S15SFCDLSM	86 - 96cm	66 - 74cm	80cm	78cm
M	S15FMED	S15SFCMED	S15SFCDLMED	96 - 108cm	74 - 80cm	90cm	86cm
L	S15FLRG	-	S15SFCDLLRG	108 - 130cm	80 - 88cm	98cm	94cm
XL	S15FXLRG	-	S15SFCDLXLRG	130 - 144cm	88 - 100cm	114cm	103cm
XXL	S15FXLRG	-	S15SFCDLXXLRG	144 - 160cm	96 - 104cm	120cm	107cm

