


Healthcare
furniture,
enhancing
lives.

Caremed Alrick.

M9000 Mental Health
Mattress

4-Way Turn, Fully Sealed &
Enclosed, with Heavy Duty
PVC Cover



Single Core Foam Mattress, allows 4 way turn which is fully sealed and enclosed. The mattress cover is a heavy duty PVC which is a double laminated surface, reinforced with additional heavy duty stitching at stress points.

The M9000 mattress is designed to minimize the damage that may be caused by abuse in high security environments. This mattress is available in both single and king single widths.

Specifications

Overall Dimensions	Single - 1960 x 870 x 150mm King Single - 1960 x 1050 x 150mm
Product codes	M9000 (Single Width) M9000KS (King Single Width)
SWL	150kg
Warranty	2 yr

FEATURES

- Single Core Foam - 4 Way Turn
- Heavy Duty PVC Cover
 - Exceptional tear resistance
 - Hypoallergenic/non-irritating
 - Anti-Static
 - Fluid proof
 - Fire retardant
 - Anti-microbial protection
 - Stain resistant
 - Fire retardant
 - Easy to maintain
 - Meets and exceeds infection control requirements
- Fully Sealed & Enclosed Cover
 - Ultimate core durability and protection, protects foam core against willful damage
 - 100% waterproof seams to prevent ingress of fluid and infection
 - High frequency sealed and welded seams
- Suitable for;
 - High Security Environments
 - Correctional Centres
 - Institutional Facilities
 - Evacuation Mattress
 - Hostels & Dormitories

