

MultiFit

The Oxford® MultiFit sling is a general-purpose sling with deluxe 'paddle' style legs for increased comfort and support. Integrated low-friction rip-stop nylon to the leg and spine sections ensures easy application and helps to reduce the risks associated with tissue viability.

- ✓ Optional detachable head support
- √ Colour coded binding indicates sling size
- Positioning handles allow easy patient repositioning
- ✓ Coloured attachment loops promote correct/safe positioning
- √ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs


















Chair to Chair



Bed to Chair/Chair to Bed



Chair to Commode/Toilet



Transfer to Bath



Floor to Bed/Chair (Head Support recommended)



Rehab/Assisted Standing

Size		Part No.	SWL
1	XS	SL1580	227 kg
2		SL1581	227 kg
3	S	SL1582	227 kg
4		SL1583	227 kg
5	M	SL1584	227 kg
6		SL1585	227 kg
7		SL1586	227 kg
8		SL1587	227 kg
9	L	SL1588	227 kg
10		SL1589	227 kg
11	XL	SL1590	284 kg
12		SL1591	284 kg







