

Perfect Positioning, Ultimate Comfort

RANGE OVERVIEW





About Finch Healthcare

Finch Healthcare was established to help improve Assistive Technology equipment outcomes in Australasia. Born of a quest for clinical excellence, our team brings equipment with proven advantages and advancements to our carefully chosen dealers. Based in Port Melbourne, VIC, we have great distribution options along with warehousing for ample stock.



About Cura Seating

Cura Seating is a custom seating manufacturer based in Peterborough UK. Their current range of products tell a story of over 40 years of experience in designing and manufacturing Care Chairs for challenging medical conditions, and with their motto "anything's possible", new models and features are frequently launched. Cura Seating also works closely with healthcare professionals in their home market which ensures all product development is results-driven.

Cura Legacy Chair



The Cura Legacy is a cost effective, highly adaptive Care Chair that has Tilt in Space, Backrest Recline and Elevating Leg rest as standard. Optional upgrades to postural support, footplates, tray and electric option along with a great weight capacity make it an excellent choice for many neurological conditions. Coming in 3 useful sizes, the Infant Legacy especially stands out in its field as unique option for children. This range is stocked in attractive colours, but a large range of options in fabrics and colours are available on special order.

CHAIR OPTIONS:

- ✓ Electric operation
- ✓ Lap tray
- ✓ Seat cushion upgrades
- ✓ Footplate padding
- ✓ Large fabric Selection

SUITABLE FOR:


Advanced Huntington's, Parkinson's, Alzheimer's, MS, MND, Paraplegia, Osteoarthritis, Stroke, Elderly and Neurological Debilitating Conditions.



Model	Infant	Kinder	Adult
Parameter	Inch/cm	Inch/cm	Inch/cm
Seat Width Range	6/15 to 16/40.6	10/25 to 16/40.6	16/40.6 to -22/56
Seat Depth Range	8/20.3 to 15/38.0	12/30.5 to 21/53	16/40.6 to -22/56
Backrest Height	19/48.26	26/66	30/76.2
Seat to Armrest Height	3/7.6 to 14/35.6	3/7.6 to 14/35.6	7/18 to 10/25
Floor to Seat Height	22/56 to 25/64	21/53 to 25/64	16/40.6 to -22/56
Seat to Footplate	8/20.3 to 25/63.5	8/20.3 to 25/63.6	8/20.3 to 25/63.6
Maximum User Weight	127kg*	127kg*	127kg*

* Custom dimensions and weight capacities available.

Cura Contour Chair



The Cura Contour chair is unashamedly the most adjustable positioning chair available. With its Oyster backrest system that boasts 25 individual adjustment points, even extreme postural deformities are readily accommodated without further customisation! Coming in 2 sizes, this chair has ample allowance for growth, making it a very viable investment. This range is stocked in attractive colours, but a large range of options in fabrics and colours are available on special order.

CHAIR OPTIONS:

- ✓ Electric operation
- ✓ Lap tray
- ✓ Seat cushion upgrades
- ✓ Footplate padding
- ✓ Large fabric selection

SUITABLE FOR:

MS, MND, Paraplegia, Scoliosis, CP, Lordosis, Kyphosis, Elderly, Neurological Debilitating Conditions.

Model	Kinder	Adult
Parameter	Inch/cm	Inch/cm
Seat Width Range	12/30.5 to 18/45.7	16/40.6 to 22/55
Seat Depth Range	12/30.5 to 20/50.8	16/40.6 to 22/56
Backrest Height	26/66	30/76.2
Seat to Footplate	8/20.3 to 25/63.6	8/20.3 to 25/63.6
Seat to Armrest Height	4/10 to 6/15	7/18 to 10/25
Floor to Seat Height	19/48.3 to 22/56	16/40.6 to 22/56
Maximum User Weight	127kg*	127kg*

* Custom dimensions and weight capacities available.



Cura Multi Porta Care Chair



The Cura Multi Porta Riser-Lift Chair has been developed with progressive conditions in mind, incorporating the following individually operating functions: Tilt in Space, Backrest Recline and Elevating Leg rest, along with a feature that enables the seat to rise vertically (seat elevation). The postural support in the seating, and footplate options are also readily upgraded. This model is stocked in attractive colours, but a large range of options in fabrics and colours are available on special order.

CHAIR OPTIONS:

- ✓ Lap tray
- ✓ Seat cushion upgrades
- ✓ Footplate padding
- ✓ Large fabric selection

SUITABLE FOR:

Palliative Care, Parkinson's, Alzheimer's, MS, MND, Osteoarthritis, Stroke, Elderly, Neurological Debilitating Conditions.



Model	Adult
Parameter	Inch/cm
Seat Width Range	16/40.6 to 24/61
Seat Depth Range	17/43.2 to 22/56
Backrest Height	30/76.2
Floor to seat height Range	17/43.2 to 25/64
Seat to Armrest Height	3/7.6 to 14/35.6
Maximum User Weight	127kg*

* Custom dimensions and weight capacities available.

Cura Multi Millennium Care Chair



The Cura Multi Millennium Riser-Lift Chair has been developed for progressive conditions, and is intended as a stationary chair. The following individually operating functions come as standard: Tilt in Space, Backrest Recline, Elevating Leg rest along with a feature that enables the seat to rise vertically (seat elevation). The postural support in the seating, and footplate options are also readily upgraded. This model is stocked in attractive colours, but a large range of options in fabrics and colours are available on special order.

CHAIR OPTIONS:

- ✓ Lap tray
- ✓ Seat cushion upgrades
- ✓ Footplate padding
- ✓ Large fabric selection

SUITABLE FOR:

Parkinson's, Alzheimer's, MS, MND, Osteoarthritis, Stroke, Elderly, Neurological Debilitating Conditions

Model	Adult
Parameter	Inch/cm
Seat Width	16/40.6 to 24/61
Seat Depth Range	17/43.2 to 21/53.3
Backrest Height	30/76.2
Floor to seat height Range	13/33 to 25/64
Seat to Armrest Height	3/7.6 to 14/35.6
Maximum User Weight	127kg*

* Custom dimensions and weight capacities available.



Innovation to Enable



Sales Enquiries

Call 1300 251 545

Email sales@finchhealthcare.com.au

Visit finchhealthcare.au

45 Bayside Avenue, Port Melbourne

Open Hours: Mon-Fri 8am – 5pm

SHOWROOM VISITS BY APPOINTMENT ONLY

SCAN ME

