

# Metzeler Mattress

## Frequently Asked Questions

### **What sizes are available?**

The mattresses are available in single, king single, double, queen, & split queen, and in soft or medium firmness.

### **Who is the mattress best suited to?**

The Metzeler mattress offers optimal comfort and support for individuals with a disability, including those with Central Nervous System (CNS) disorders, post-surgery or cancer treatment, and conditions that may impact thermoregulation such as CP, MND, and MS. Moreover, it provides exceptional comfort and quality sleep to those who haven't had success with standard static foam or alternating air mattresses.

### **How does the mattress maintain or support mobility and function?**

The reinforced edges make it easier to get in and out of the bed. The vertical tubes also provide stability when repositioning the body.

### **What makes the Metzeler mattress different to other medical foam mattresses?**

It is the only static mattress to offer vertical and horizontal pressure support. This is achieved through an outer layer of castellated foam which assists air flow and reduces positive resistance. It also features an inner core with two different tube types (firmness). This provides even pressure distribution, supporting the natural curves of the spine.

### **What is meant by multizone system and what are its benefits?**

The mattress contains 7 zones of differing firmness to correctly support the body and provide optimal pressure redistribution.

### **Given the mattress is foam, is it hot? Does it increase sweating?**

Metzeler's patented Rubex tube system optimises ventilation, breathability and enhances thermoregulation. Heat and moisture are transferred through and out of the mattress providing a healthy microclimate between the mattress surface and client's body.

### **Is the mattress environmentally friendly?**

The Metzeler mattress is made from natural non-toxic environmentally friendly raw materials. The main component is sunflower oil.

### **Does the mattress have a warranty?**

All Metzeler mattresses come with an eight year warranty.



**Call: 1300 251 545**

Visit: [finchhealthcare.au](http://finchhealthcare.au)

Email: [sales@finchhealthcare.au](mailto:sales@finchhealthcare.au)

45 Bayside Ave, Port Melbourne 3207