

Healthcare
furniture,
enhancing
lives.

2300T

2300 Tilt Series Bed



**Caremed
Alrick.**



AUSTRALIAN MADE
AND OWNED

23T-R2

2300T



The 2300T Series bed is the most versatile Hi-Low and tilt bed available for the healthcare industry. Featuring up to 300kg SWL sturdy frame, therefore being suitable for patients with a range of weight requirements, this bed offers sizes from single up to king widths. The 2300T bed is excellent for in-home, and facility care as it has many options available making this bed aesthetically pleasing and practical for both the patient and carers.

FEATURES

- Electric Hi-Low
- Manual Vascular Support
- Back lit Handset
- APO Grey & Ironstone Powder Coated Frame
- Curved Natural Oak Head & Footboards
- Electric Kneebreak
- Electric Backrest
- Forward and Reverse Trendelenburg
- Bluetooth Connectivity

****handset must be ordered separately****



Pull Out Extension



Auto Chair Position



125mm Braking Castors



Underbed Light

SPECIFICATIONS

2300 Tilt	Single	King Single	Extra Wide	Double	Queen	King
Overall Dimensions (mm)	2040 x 920	2040 x 1090	2040 x 1220	2040 x 1390	2040 x 1520	2040 x 1820
Sleeping Surface (mm)	2000 x 900	2000 x 1070	2000 x 1200	2000 x 1370	2000 x 1500	2000 x 1800
SWL (kg)	300	300	300	300	300	300
High/Low Position (mm)	830/390	830/390	830/390	830/390	830/390	830/390
Back Angle (degrees)	77	77	77	77	77	77
Knee Angle (degrees)	40	40	40	40	40	40
Reverse Angle (degrees)	16	16	16	16	16	16
Trendelenburg Angle (degrees)	16	16	16	16	16	16
Product Codes	23ST-000100	23KST-000100	23EWT-000100	23DT-000100	23QT-000100	23KT-000100

WARRANTY: Frame - 10 years, Electrics - 5 years & Castors 1 year.

2300T

OPTIONAL FEATURES

Added at time of manufacture

- Custom Powder Coat Colours (as shown)
- Head & Footboard Options
- Attendant Control
- Electric Vascular Support



Bluetooth Handset



Battery Back Up



CPR Backrest



Central Locking Castors



OPTIONAL HEAD & FOOTBOARDS

Custom options available

Laminate Curved Style



● Natural Oak

Upholstered Jersey Headboard & Plain Footboard



● Ash
● Ocean

Timber Eltham Style



● Burnt Walnut
● Beech
● Saddle



ACCESSORIES

- Mattresses
- Cot Pro
- T-Bar Buffer



400HR
Horizontal Side
Rails



380VR
Vertical Side
Rails



CLP400
Clamp on
Support Rail



CGH
Combo Grab
Handle



SE-SH
StandEzy - Straight



SE-HH
StandEzy - Hook
Handle



ASHP/H
Self Help Pole with
Triangle



AIVP
IV Pole


2300T



Can't find what you're looking for? Contact us to discuss your custom options.



Pick your colour



Choose your width



Measure your required length



Create your style



We build for legends



Connect with us

Distributed by:

