

Healthcare
furniture,
enhancing
lives.

Caremed Alrick.

NINO Series Bed

The Nino Series allows children to perceive the world around them while at the same time providing protection and security. The transparent Perspex merges with the room allowing occupants to have an unrestricted view of their immediate surroundings yet still protects against falls and injuries. Also available in different powdercoat colours, or with higher door elements or doors on both sides.



** Colours of the Nino in this image are a custom design these are not the standard powdercoat and endboard colours **

Configuration Possibilities



Options Available



Doors on both sides (only available on standard height doors)

Higher doors

Central locking

Features



Standard features also include:

- Doors one side (750mm high perspex sides)
- 125mm braking and directional castors
- Battery backup
- Underbed lighting
- Backlit handset
- White laminate

Specifications

| | | |
|-----------------------------|-------------------------------|-------|
| Overall dimensions | 2200 x 1014 mm | |
| Sleeping surface dimensions | 2000 x 870 mm | |
| SWL | 200kg | |
| NET Bed weight | 140kg | |
| Low Position | 390mm | |
| High Position | 830mm | |
| Underbed clearance | 100mm (at lowest point) | |
| Angles (°) | Back | 70 |
| | Knee | 44 |
| Warranty | Frame | 10 yr |
| | Electric | 5 yr |
| | Castors | 1 yr |
| | Handset | 1 yr |
| Optional features | Powdercoat colours | |
| | Higher doors (1100mm) | |
| | Doors both sides (750mm only) | |
| | Central locking castors | |
| Product codes | NINO | |

** LOAD CAPACITY HAS BEEN TESTED IN ACCORDANCE WITH AS/NZS 3200.2.38:2007 **

Powdercoat Colours Available

- Please ensure you state your 2 preferred colours for the Nino. See below for details.
- Pricing includes 2 Powdercoat colours from the below selections.
- If there is no colours stated on the PO, Nino will be supplied in our standard APO Grey / Ironstone tones.

