Healthcare furniture, enhancing lives.



Plexa Series Bed

















The Plexa Series Bed is designed to provide safe and comfortable sleeping arrangements for individuals with complex needs. This bed includes full-length rails to provide security, featuring a locking system, ensuring the rails stay in place while the user is sleeping.

An additional feature is the height adjustment making transfers in and out of bed with ease, designed with a durable frame and available in various sizes to accommodate the user's needs.

FEATURES

- APO Grey Powdercoated Frame
- Auto-regression Back Rest
- Trendelenburg
- 75mm Castors
- Vertical Lift
- Extremely Low Minimum Position







Bluetooth connectivity Timber Drop ***handset must Down Rails be ordered separately***



Full Length



Chair Position at the push of a button

SPECIFICATIONS

PLEXA	Single	King Single	Double
Overall dimensions (mm)	2150 x 1040	2150 x 1210	2150 x 1510
Sleeping surface (mm)	2000 x 920	2000 1070	2000 x 1370
SWL (kg)	250	250	250
Low/high position (mm)	170/725	170/725	170/725
Back angle (degrees)	62	62	62
Knee angle (degrees)	35	35	35
Reverse angle (degrees)	15	15	15
Trendelenburg angle (degrees)	15	15	15
Product Codes	PLEXA-S-00000	PLEXA-KS-00000	PLEXA-D-00000

WARRANTY: Frame - 10 years, Electrics - 5 years & Castors -1 year.



T: 1800 282 001

E: sales@caremed.healthcare W: www.caremed.healthcare





OPTIONAL FEATURES

Added at time of manufacture

- **Dual Locking Castors** 0
- **Underbed Light** 0
- Battery Back Up



Bluetooth Handset



Custom Powdercoat Colours



Head & Footboard Styles (as shown in Jersery)



Head & Footboard Colour (as shown in Stone)



OPTIONAL HEAD & FOOTBOARDS

Custom options available

Laminate Square Style Head & Footboard



Natural Oak

Upholstered Plain Padded Head & Footboard



Ocean

ACCESSORIES

- **Emergency Power Supply**
- **ASHP**
- Mattress





Reading Light



IV Pole



EN9-TBR



Cot Side **Protectors** with Mesh Insert



T: 1800 282 001

Alrick. E: sales@caremed.healthcare W: www.caremed.healthcare







Can't find what you're looking for? Contact us to discuss your custom options.



colour



Choose your width



Measure required length



Create your style




We build for legends



Connect with us

Distributed by:





T: 1800 282 001 Alrick. E: sales@caremed.healthcare W: www.caremed.healthcare

