

Healthcare
furniture,
enhancing
lives.

SECURA 2

Secura 2 Mental Health Bed



SECURA 2




The Secura 2 provides dignity and comfort for all mental health patients, with its impressive safe working load of 300kg, designed specifically to reduce stress on carers by relieving physical strain. Secura 2 is fitted with heavy-duty actuators and concealed wiring to ensure reliable and efficient functioning. Suited for low to medium care patients who may have clinical care needs with reduced anti-ligature features.

FEATURES

- Electric Hi-Low
- Enclosed anti-ligature area
- All bolts secured with Loctite
- Ironstone Powdercoat



Embedded footboard control



Concealed wiring



Braided electrical power supply



Anti-tamper screws



Secure bolt on head & footboards



Curved Natural Oak head & footboard



Enclosed mattress supports



125mm braking castors

SPECIFICATIONS

Secura 2	Single	King Single
Overall dimensions (mm)	2070 x 920	2070 x 1090
Sleeping surface (mm)	2000 x 900	2000 x 1070
SWL (kg)	300	300
Low/high position (mm)	390/830	390/830
Product Codes	SEC2-S-000010	SEC2-KS-000010

WARRANTY: Frame - 10 years, Electrics - 5 years & Castors 1 year.

SECURA 2

OPTIONAL FEATURES

Added at time of manufacture

- Knee break
- Trendelenberg
- Bluetooth handset
- Custom powdercoat colours
- Extra length



Back rest
(as shown)



Battery
back up



Central
locking
castors



Key power
lockout
(as shown)



OPTIONAL HEAD & FOOTBOARDS

Custom options available

Laminate

Curved Style

(25mm headboard & 32mm footboard)



Natural Oak



Oyster Grey

ACCESSORIES

Mental Health Mattress

M9000MH (single)

M9000KSMH (king single)




SECURA 2



Can't find what you're looking for? Contact us to discuss your custom options.



Pick your colour



Choose your width



Measure required length



Create your style



We build for legends



Connect with us

Distributed by:

