

Healthcare  
furniture,  
enhancing  
lives.

# SECURA 3

Patented Secura 3 Mental Health Bed




\*\*\* Pictured with optional back rest \*\*\*

# SECURA 3




\*\*\* Pictured with optional back rest \*\*\*



The Secura 3 is the ultimate safe and secure solution for your facility, eliminating access to the Patented bed mechanisms, with the bed surround maintaining the concealment in both, Hi and Low positions, the lower bed section sitting stationary with the bed. The sides of the bed are recessed to allow for nurses to kneel on the bed without the obstruction of a bed surround in the shin area. Secura 3 is the solution to your WHS concerns for the carers and cleaners.

## FEATURES

- Concealed castors



Embedded footboard control




Brake lever access



Natural Oak, 13mm Compact Laminate



Ironstone powdercoat



Lowered section for kneeling on bed



Head & footboard, lower than mattress height



Braided electrical power supply (700mm in length)



Anti-tamper screws

## SPECIFICATIONS

Secura 3	Single	King Single
Overall dimensions (mm)	2030 x 930	2030 x 1100
Sleeping surface (mm)	2000 x 900	2000 x 1070
SWL (kg)	250	250
Low/high position (mm)	400/700	400/700
Product Codes	SEC3-S-010010	SEC3-KS-010010

WARRANTY: Frame - 10 years, Electrics - 5 years & Castors 1 year.

# SECURA 3

## OPTIONAL FEATURES

Added at time of manufacture

- Kneebreak
- Custom powdercoat colours
- Longer sleeping surface
- Custom length power cable




Back rest  
(as shown)



Battery  
back up



Bluetooth  
Controls



Full height  
head &  
footboard



## OPTIONAL SURROUND

Custom options available

### Compact Laminate

(13mm Compact Laminate)



Natural Oak



Oyster Grey

\*\*\* Pictured with optional back rest \*\*\*

## ACCESSORIES

### Mental Health Mattress

M9000MH (single)

M9000KSMH (king single)



# SECURA 3



\*\*\* Pictured with optional back rest \*\*\*

Can't find what your looking for? Contact us to discuss your custom options.



Pick your colour



Choose your width



Measure required length



Create your style



We build for legends



Connect with us

Distributed by:

