


# Good Night, Better Day



**GET THE PERFECT SLEEP**  
Learn more inside



# Get the perfect sleep with Metzeler Mattress

We spend more than one third of our life sleeping, moving into up to 60 different positions per night. Therefore it's important you choose a mattress that delivers you a quality sleep to ensure you wake up feeling refreshed.















The Metzeler mattress will give you with just that - a perfect night's sleep. Incorporating patented Rubex™ Duo Cell technology, the mattress offers superior comfort and support which provides optimal pressure distribution, ventilation and spinal alignment.



## Featuring Rubex Duo Cell Patented Technology

Rubex Duo Cell consists of approximately 50% small and large pores. The large pores reduce pressure and the small ones provide soft pressure support. It also provides excellent air permeability and moisture transport to help create a constant bed temperature. The foam reacts only to body pressure providing constant support over a wide temperature range (-5°C - +40°C.)

## Metzeler Product Range


Mattress Size	Firmness	Dimensions	Product Code
  <b>Single</b>	Medium	900x2000x180mm	729595-ADM090
  <b>Single</b>	Soft	900x2000x180mm	732758-ADS090
  <b>King Single</b>	Medium	1050x2000x180mm	729611-ADM105
  <b>King Single</b>	Soft	1050x2000x180mm	732757-ADS105
  <b>Double</b>	Medium	1350x2000x180mm	721972-ADM135
  <b>Queen</b>	Medium	1520x2000x180mm	721970-ADM152
  <b>Split Queen</b>	Medium	760x2000x180mm	721967-ADM076





### Waterproof & breathable cover

The waterproof and breathable visco-elastic cover assists with moisture management, while being easily removable and washable up to 95°C.



### Optimal temperature management

The tubular structure and a choice of soft or medium foam castellated overlay on top and bottom of the inner core provide optimal horizontal and vertical cross-ventilation thermo-regulation needs.



### Durable reinforced mattress edge

Reinforced foam edges assist carers and physios with patient transfer, while minimizing the potential for rolling back into or off the mattress.




### Multizone support system

The core of this mattress has a unique tubular structure that is zoned to support the spine and provides optimal pressure redistribution.



### Environmentally friendly design

Unlike most foam mattresses made from petroleum byproducts, Metzeler mattresses use natural, non-toxic materials like sunflower oil and other environmentally-friendly alternatives.





## About Finch Healthcare



### Quality Equipment

We supply retail partners with quality equipment for individuals with special needs, disabilities or the elderly.




### Industry Knowledge

Our team is knowledgeable about disability care and committed to outstanding customer service.



### Retailer Support

We offer expert training and support to ensure the end client receives the best possible equipment solution.



**Call: 1300 369 096**

Visit: [finchhealthcare.au](https://finchhealthcare.au)

Email: [sales@finchhealthcare.au](mailto:sales@finchhealthcare.au)

45 Bayside Ave, Port Melbourne 3207