

CATALYST A & AC





INNOVATIVE AND PATENTED FRAME DESIGN

Patented frame concept uses intelligent geometry to create a significantly stronger, lighter and better performing chair.



OFFSET AXLE PLATE

The Catalyst 4 series vertically adjustable axle plate is mounted forward of the back frame. This small change can significantly improve efficiency.



SWING AWAY FOOTRESTS

Enjoy the simplicity and strength of Catalyst's patented swing-away footrest hangers. The tapered engagement keeps the footrests snug and rattle free. Comes with a lifetime warranty.

EASY TO TRANSPORT. EASY TO PUSH.

The Catalyst 4 & 4C are the lightest, most durable and easiest to use K4 folders on the market.



SPECIFICATIONS

Transport Weight:	16 lbs.*	7.25 kg*
Seat Width:	14-22 in.	36-56 cm
Seat Depth:	14-20 in.	36-51 cm
Front Seat Height:	13.5-21 in.	34-53 cm
4 Rear Seat Height:	12.5-19 in.	32-48 cm
4C Rear Seat Height:	12.5-20 in.	32-51 cm
Top Frame Tubing Size:	1 in.	2.54 cm
Bottom Frame Tubing Size:	1.125 in.	2.86 cm
Weight Capacity:	250/350 lbs.	113/158 kg

* 16"x16" Frame, 5"x1" Poly Casters, Height Adjustable T-Arm Receivers. Footrest and Rear Wheels removed.