

# CATALYST X



# 5VX

# IGNITING THE EVOLUTION

ENDLESS **STYLE**. EXPECTED **QUALITY**. ESSENTIAL **VALUE**.



From the logo to the frame, Catalyst 5VX has been redesigned with a modern, minimalist approach. Enhanced designs let you optimize your ride.

**THE RIDE YOU'VE COME TO EXPECT FROM CATALYST.**

## **BETTER STYLE**

Be bold. Be you. With thousands of combinations between frame color, anodizing package, options and accessories, your style begins with you.



## **BETTER ADJUSTMENT**

Center of Gravity adjustment and ride performance go hand in hand. The redesigned axle plate and rear frame give you more adjustment over a greater range. Find your fit, ease your propulsion and optimize the ride you've come to expect with Catalyst.

Adjustable Camber gives you up to 4° of camber in 1" increments. A simple adjustment and convenient option you can optimize now or later without any additional parts.



# MORE OPTIONS & ACCESSORIES

## SWING-IN ANTI-TIPS

The new swing in anti-tips make it easy for the user to operate. By simply pushing down on the ergonomically designed handle you can swing them in, out of the way.



## TRANSIT WHEELS

Integrated into the rear frame, Transit Wheels allow you to remove your rear wheels for easier access through narrow spaces.



## FENDER SIDE GUARDS

Choose between carbon fiber and foam padded aluminum side guards with fender to protect your clothing from the road while adding a sleek look to your chair.



### SPECS

Transport Weight*	13.9 lbs.*	6.3 kg
Seat Width	14-22"	36-56 cm
Seat Depth	14-20"	36-51 cm
Front Seat Height	13.5-21"	34-53 cm
Rear Seat Height	11.5-20"	29-51 cm
Top Frame Tubing Size	1"	2.54 cm
Bottom Frame Tubing Size	1.125"	2.86 cm
Weight Capacity	300 lbs.	136.08 kg
Ansi Resna WC-19 Transit Approved	Yes	Yes

\* 16"x16" Frame, 5"x1" Poly Casters, Height Adjustable T-Arm Receivers. Footrest and Rear Wheels removed.



Stevens Point, WI 54481 | 715-254-0991 | [www.kimobility.com](http://www.kimobility.com)

DCN0319.0