

Little Wave





CRESCENT CG SYSTEM

Hidden inside your Little Wave, the simple, yet ingenious, Crescent CG system keeps the side frame free of excess hardware while providing infinite adjustability. It is this infinite adjustability that allows for CG positioning anywhere along the side frame without interference with the backrest hardware or aftermarket seating system.

MULTI-ADJUSTABLE FLIP-UP FOOTPLATE

This multi-adjustable footrest system allows for an infinite number of positions so it can be repositioned as the rider grows. The aluminum footplate with carefully radiused edges, will withstand the wear and tear of any army...or a 5 year old! In addition, the easy access strap allows the rider to lift up the footplate without bending down.

GROWING WITH YOU

The Little Wave grows in depth and also in width by using frame Xpanders that maintain the frame's structural integrity far better than bolt-on cross bars. Knowing you will need to grow your Little Wave XP, this version includes one free growth kit consisting of everything you need to grow the frame!

Imagine... a wheelchair that has the ability to move easily within its environment; where the wheels are positioned perfectly so young people can use all of their strength and leverage to their advantage.

Imagine... a wheelchair that has the ability to change and grow along with the person using it yet remains extremely light, simple to use and easy to transport.

It's not imaginary. It's reality.
Little Wave.



SPECIFICATIONS

Transport Weight:	12.5 lbs.*	5.67 kg*
Seat Width:	10-16 in.	25-41 cm
Seat Depth:	10-16 in.	25-41 cm
Front Seat Height:	11.5-20.5 in.	29-52 cm
Rear Seat Height:	11.5-20 in.	29-51 cm
Tubing Size:	1 3/8 in.	3.49 cm
Weight Capacity:	165 lbs.	75 kg
Ansi Resna WC-19 Transit:	Yes	Yes

* 12"x12" Frame with Adjustable Height Back - 12", 4" Poly Casters, Angle Adj. Footplate and Standard Upholstery. Armrest and Rear Wheels Removed.