



Raz-AP

Attendant Propel

Rehab S hower Commode

The **Raz-AP** (Attendant-Propel) rehab s hower commode chair comes standard with 4" (10 cm) of height adjustment, tension-adjustable back upholstery.



FEATURES

Extensive range of options available

Flip-up padded arms

Exclusive to Raz, Ischial Pelvic Alignment System (IPAS), provides 2" (5 cm) of fore/aft or 1" (2.5 cm) of lateral seat adjustment

Stainless steel frame with 4" (10 cm) of height adjustment



Tension-adjustable back upholstery

Adjustable back angle

Dual-locking (swivel and roll) casters

ALSO AVAILABLE

- Long frame - 3" (8 cm) longer for tall clients
- Low frame - 4" lower seat-to-floor height
- 16" and 20" (40 and 51 cm) widths
- Recliner frame

SPECIFICATIONS	Raz-AP (Attendant Propelled)
WEIGHT CAPACITY	350 lbs (158 kg)
FRAME WIDTH (seat and back frame)	18" (46 cm)
OVERALL WIDTH (with standard arm pad)	22.5" (57 cm)
OVERALL LENGTH (with footrest)	33.75" (86 cm)
SEAT WIDTH (commode seat)	18" (46 cm)
SEAT DEPTH (front of back post to front of seat)	16" or 19" (40 or 48 cm)
CLEARANCE HEIGHT UNDER SEAT	15.5" - 19.5" (39 - 50 cm) with pan brackets 16.25" - 20.25" (41 - 51 cm) without pan bracket
ARMREST HEIGHT ABOVE SEAT (17" Adjusta-Back)	6.5" - 9.5" in ½" increments (17 - 24 cm)
SEAT-TO-FLOOR HEIGHT	20" - 24" in 1" increments (51 - 61 cm)
PRODUCT WEIGHT (with footrests)	38 lbs (17.25 kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.

MANUFACTURED BY:

RAZ Design Inc. 22 Howden Road, Toronto, Ontario Canada,

M1R 3C7 T: (877) 720-5678 E: info@razdesigninc.com

RAZSAIG002r01 Rev. Date: 2017-02-27

Raz Design products are distributed in Australia by:
Paragon Mobility

3/10 Tapnor Crescent, Brendale, QLD, 4500 Australia
Tel: 13000 652 382
Fax: (07) 3112 6125
www.paragonmobility.com.au

