

Find Your Freedom



RANGE OVERVIEW

ARTG ID 429750



Mountain Trike

The Mountain Trike is a manual, all-terrain wheelchair with a unique lever drive system for clean hands. It's versatile for daily outdoor use, adaptable for children, and can convert to MT Push or eTrike. Comes with a 3-year warranty.

Feature	Description
Frame	Aluminium Frame
Suspension	Rock Shox monarch R Air shocks
Brakes	All weather SRAM Level hydraulic disc
Footrest	Adjustable with foot straps
Support	Moulded cushion, seatbelt (optional)
Length	132 cm
Seat Width	17 in
Seat Angle	20 degrees
Seat Cushion	Removable foam cushion
Seat Height	46.5 cm (from ground)
Seat Length	44.7cm - 48.7 cm
Footrest Height	41cm - 48 cm
Wheel Size	24 in
Overall Width	74 cm
Weight	20 kg
Back Height	37 cm (adjustable)
Rider Max	100 kg, 6'4"

Mountain Trike Push

The MT Push is an attendant all-terrain wheelchair with adjustable push handles, ideal for inclusive access. Developed with the National Trust, it's perfect for country walks, city breaks, and festivals. Convertible to lever drive Mountain Trike. Comes with a 3-year warranty.

Feature	Description
Frame	Aluminium Frame
Suspension	Rock Shox monarch R Air shocks
Brakes	All weather SRAM Level hydraulic disc
Footrest	Adjustable with foot straps
Support	Moulded cushion, seatbelt (optional)
Length	132 cm
Seat Width	17 in
Seat Angle	20 degrees
Seat Cushion	Removable foam cushion
Seat Height	46.5cm
Seat Length	44.7 - 48.7 cm
Footrest Height	41 - 48 cm
Wheel Size	24 in
Overall Width	74 cm (78cm push rims)
Weight	20 kg
Back Height	42 cm (adjustable)
Rider Max	100 kg, 6'4"



Mountain E-Trike

The eTrike is an electric-assist all-terrain wheelchair, perfect for those needing assistance or wanting to tackle challenging terrain. It offers manual, power assist, and hybrid modes. Range is 8-20 miles, depending on conditions. Comes with a 3-year warranty.

Feature	Description
Frame	Aluminium Frame
Suspension	Rock Shox monarch R Air shocks
Brakes	All weather SRAM Level hydraulic disc
Footrest	Adjustable with foot straps
Support	Moulded cushion, seatbelt (optional)
Length	132 cm
Seat Width	17 in
Seat Angle	20 degrees
Seat Cushion	Removable foam cushion
Seat Height	46.5cm
Seat Length	44.7 - 48.7 cm
Footrest Height	41 - 48 cm
Wheel Size	24 in
Overall Width	74 cm
Weight	25 kg inc battery
Back Height	42 cm (adjustable)
Rider Max	100 kg, 6'4"
Motor	250W, 36V Hub Motor
Battery	12.8Ah, 36V Lithium Ion
Speed	max 8mph (12.8kmh)



Mountain E-Trike Push

The ePush is an attendant all-terrain wheelchair with electric assist for easier pushing over tough terrain. It can be manually pushed or switched to electric mode. Ideal for inclusive outdoor adventures. Comes with a 3-year warranty, battery with 12 months.

Feature	Description
Frame	Aluminium Frame
Suspension	Rock Shox monarch R Air shocks
Brakes	All weather SRAM Level hydraulic disc
Footrest	Adjustable with foot straps
Support	Moulded cushion, seatbelt (optional)
Length	132 cm
Seat Width	17 in
Seat Angle	20 degrees
Seat Cushion	Removable foam cushion
Seat Height	46.5cm
Seat Length	44.7 - 48.7 cm
Footrest Height	41 - 48 cm
Wheel Size	24 in
Overall Width	74 cm (78cm push rims)
Weight	25 kg inc battery
Back Height	42 cm (adjustable)
Rider Max	100 kg, 6'4"
Motor	250W, 36V Hub Motor
Battery	12.8Ah, 36V Lithium Ion
Speed	max 3.5mph (5.6kmh)

Choose Your Perfect Mountain Trike

Features	eTrike	Mountain Trike	MT Push	ePush
Direct Steering	✓	✓	✗	✗
Air Suspension	✓	✓	✓	✓
Disc Brakes	✓	✓	✓	✓
Drive Levers	✓	✓	✗	✗
Self Propelled	✓	✓	✗	✗
Attendant Propelled	✗	✗	✓	✓
Configurable Controls	✓	✓	✓	✓
Electric Assist	✓	✗	✗	✓
Electric Assist	✓	✗	✗	✓

ARTG ID 429750

About Mountain Trike

Based in the United Kingdom, Mountain Trike is a pioneering business dedicated to revolutionising off-road mobility solutions. Specialising in all-terrain wheelchairs, they combine cutting-edge technology with robust design to provide unparalleled freedom and independence for users.

Committed to quality and innovation, Mountain Trike integrates features like hydraulic disc brakes, air suspension, and ergonomic seating to ensure safety and comfort. Their mission is to empower individuals with mobility challenges, enabling them to explore the great outdoors without limitations.

Customise Your Trike

Your Mountain Trike can be customised with an extensive range of accessories to allow you to explore the trails. Chat to your Mountain Trike dealer about how you can get the best fit and experience.



Endless Colour Options

Mountain Trike is available in a range of standard ready to ship colours, but offers an extensive range of custom colours on request.



Sales Enquiries

Call 1300 251 545
sales@finchhealthcare.com.au
 Visit finchhealthcare.au



Call: 1300 369 096
 Visit: urgoform.com.au
 Email: sales@urgoform.com.au
 Office: 5-7 Tarkin Ct, Bell Park VIC 3215

