

permobil

Elevate play

M300  
PS JR.





Highly adjustable seat that grows as fast as you do.



### Catlink Suspension

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



### Highly Configurable


Allows you to maintain proper positioning as you grow.




### Multi-motion Armrest

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

## Shroud Colour Options



## Additional Hub Colour Options



## Specifications

Maximum User Weight	75 kg
Maximum Speed	10 km/h
Range	20 km*
Base Width	24"
Base Length	36.5"
Minimum Turning Radius	20"
Ground Clearance	3"
Recommended Battery Type	Group 24, Group 34
Weight Including Batteries	142 kg G34 / 151 kg G24
Drive Electronics	120A R-Net
Maximum Slope	9° up, 9° down <sup>1</sup>
Seat-to-Floor Height	17" (elevator) 17", 17.75", 18.5" (seat tube)

Optional Seat Elevator	8"
Tilt Options	0° – 45° (power) 0° – 45° (manual)
Manual Recline	85° – 120°
Legrest Elevation	70° (swing-away) 90° – 150° (centre mount elevating)
Backrest Height	16.5" – 22.5"
Seat Depth	10" – 18"
Seat Width	11" – 16"
Armrest Pad Lengths	10", 13"

\*Actual driving range will vary based on driving and battery conditions.  
<sup>1</sup>Refer to Owner Manual for safe driving practices.