

LECKEY

Product Range Guide

Australia



Say hello to the
Leckey family of products.



I am delighted to introduce our new Leckey product guide and celebrate our Leckey family circle.

At Leckey we believe in ‘Supporting Children’.

To everyone in Leckey, Supporting Children means working together with our family of experts – designers, therapists, parents and the broader clinical community to achieve our goals.

Our **Mission** is; To improve the quality of life of children with special needs.

Our **Vision** is; A world where each child is included and they can access their full potential.

Our **Courageous Goal** is; To prevent hip migration in children with CP by 2025.

In Leckey we see hip migration as the single biggest life event that our children may experience and has the greatest impact on their quality of life and their carers from that moment on.

Everything we do in Leckey, in partnership with our family of experts, is aimed at achieving these goals.

Our approach is simple – lead through innovation and clinical focus.

Through early intervention, childhood right up to adulthood, Leckey’s experienced team of designers and clinicians work together with our family of experts to develop exciting innovative products that make a real difference to the quality of life for the children and their families.

We aim to develop clinically and technically superior solutions, products and programs, that improve the child’s Mobility, Ability and Participation.

All of our designs undergo extensive clinical and usability trials with physiotherapists, occupational therapists, users and their families.

They are supported by in-depth clinical research studies where we work with leading universities and clinicians internationally to ensure that all Leckey products deliver on the clinical needs of the healthcare professionals and the social needs of the user.

We are committed to our Mission, Vision and Courageous Goals. We want to work with you in partnership to build our family - a community of experts with shared goals and values.

Join our family and together we will make a difference.

- The Leckey Team



Contents

Seating and Mobility 04

| | |
|------------------------|----|
| Early Activity System | 06 |
| Pal Seat | 08 |
| Squiggles Seat | 09 |
| Mygo Seat | 10 |
| Everyday Activity Seat | 12 |
| Therapy Bench | 14 |

Standing and Walking 16

| | |
|-------------------|----|
| Squiggles+ | 18 |
| Squiggles TT | 20 |
| Squiggles Stander | 22 |
| Mygo Stander | 24 |
| Horizon | 26 |
| Pronestander | 28 |
| Freestander | 29 |
| Totstander | 29 |
| MyWay | 30 |
| MyWay Pedal | 32 |

Other Products 34

| | |
|---------------------|----|
| Advanced Bath Chair | 36 |
| Shower Trolley | 37 |
| Sleepform | 38 |

Seating and Mobility

LECKEY Early Activity System

Colourful, comfortable, supportive activity system for children aged 0-48 months.



Product Features

- Offers a gentle, playful introduction to 24-hour postural management.
- Practical all-in-one kit, providing versatile support. The positioning possibilities are endless and everything can be packaged up in the duffle bag for easy transport.
- Supplied with the popular Activity Development Programme which supports parent-led therapy at home.
- Soft-touch, anti-bacterial, machine washable and tumble dry safe covers.



Key Positions

- Back Lying (Supine)
- Tummy Lying (Prone)
- Side Lying
- Floor Sitting (Long sitting)
- Hands and knees (Four point kneeling)

| Technical Details | | |
|---------------------|------------------------------|----------------------|
| One Size Only | | |
| Age | Max 36 months | |
| Weight | Max User Weight 15kg / 33lbs | |
| Stature | Max 1000mm / 39.4 inches | |
| Floor Mat | Width | "760mm / 26.3 inches |
| | Length | 950mm / 37.4 inches" |
| Sitting Support | Height | "190mm / 7.4 inches |
| | Depth | 300mm / 11.8 inches" |
| Seat width | Min | "200mm / 7.9 inches |
| | Max | 250mm / 9.8 inches" |
| Support Roll small | Width | "60mm / 2.4 inches |
| | Length | 600mm / 23.6 inches" |
| Support Roll medium | Diameter | "80mm / 3.2 inches |
| | Length | 700mm / 27.5 inches" |
| Trunk Wedge | Width | "270mm / 10.6 inches |
| | Length | 208mm / 8.2 inches |
| | Angle | 7.5 degrees" |
| Positioning Straps | Width | "80mm / 3.14 inches |
| | Length | 660mm / 25.9 inches" |



LECKEY Pal Seat

Colourful, comfortable classroom seating system providing mild support for children aged 1-12 years.

Product Features

- Two-point lap belt and side pads for lateral stability. An optional ramped cushion can be added if extra help is needed to keep the pelvis and femurs in alignment.
- Where a less upright angle is required, the backrest angle and armrests adjust to suit the child's needs, and along with the large tray or table, provide appropriate support.
- Height adjustable armrests, legs and footplate.
- Available in two bright colours in wipe clean fabric [please see website for details].
- A range of base and footrest options available, including castors, rockers and stabilisers.



Available in 2 fabric colours

| Technical Details | | | | |
|-------------------|--|---|--|--------------------------------------|
| | Size 1 | Size 2 | Size 3 | Size 4 |
| Age | 1-3 | 2-5 | 4-8 | 7-12 |
| User Height | Min 700mm [27.6"] Max 950mm [37.4"] | Min 900mm [35.4"] Max 1050mm [41.3"] | Min 1000mm [39.4"] Max 1300mm [51.2"] | Min 1200mm [47"] Max 1600mm [63"] |
| Max. User Weight | 30kg | 35kg | 40kg | 50kg |
| Seat Depth | 180mm-240mm [7"-9.4"] | 220mm-290mm [8.6"-11.4"] | 250mm-330mm [9.8"-12.9"] | 300mm-400mm [11.8"-15.7"] |
| Seat Width | 150mm-250mm [5.9"-9.8"] | 170mm-270mm [6.6"-10.6"] | 230mm-330mm [9"-12.9"] | 280mm-380mm [11"-14.9"] |
| Seat Height | 160mm-240mm [6.2"-9.4"] | 185mm-285mm [7.2"-11.2"] | 250mm-350mm [9.8"-13.7"] | 300mm-400mm [11.8"-15.7"] |

LECKEY Squiggles Seat

Suitable for use as a home chair, school chair or with a Zippie base, Squiggles meets all the postural and comfort requirements of a seating system for young children with special needs, while maintaining a fun, colourful design.

Product Features

- Designed specifically for the developmental needs of younger children- the one-piece seat base and adjustable sacral support can be used to accommodate flexion or encourage extension.
- The shape of the seat and back, the laterals with their fabric tip which wrap securely around the child and the way in which the 4-point pelvic harness integrates with the hip pads, provides a snug, comfortable fit and streamlined look.
- Being low to the ground, the Squiggles Seat promotes peer to peer interaction, making it ideal for home and nursery/school environments.
- Adjustable footplate can be tilted to the preferred angle, ensuring the feet are well supported and providing a secure base for maximum upper body function. The footplate flips away to allow for easy transfers.
- Tilt-in-space feature provides a simple way to achieve a change of position without disturbing overall posture.



Available in 3 fabric colours

| Technical Details | |
|------------------------------|----------------------------|
| One Size Only | |
| Age | 1-5 |
| User Standing Height | 750mm-1100mm [29.5"-43.3"] |
| Max. User Weight | 22kg [48.4lbs] |
| Seat Depth | 190mm-290mm [7.5"-11.4"] |
| Seat Width | 160mm-260mm [6.3"-10.2"] |
| Chest Width | 160mm-220mm [6.3"-8.75"] |
| Backrest Height | 320mm-380mm [12.6"-15"] |
| Lower Leg Length | 122mm-270mm [4.8"-10.6"] |
| Tray Height | 115mm-170mm [4.5"-6.7"] |
| Backrest Angle | Prone 10° - Recline 25° |
| Seat Height (Hi-Low Chassis) | 270mm-650mm [10"-25"] |

Squiggles Seat now able to interface with Zippie IRIS
Crash Tested Combination!



LECKEY Mygo Seat

The Mygo offers postural support whilst enabling function and mobility. It is available in three different sizes for children aged 3-16 years.



Product Features

- Second size in Leckey's family of clinically based seating systems.
- 4-point harness on the pelvic harness and pelvic cradle ensures the user stays in the correct anatomical position for optimum support.
- The flexible sacral support works with the hip guides and adjustable seat angle to provide the optimum achievable posture for the pelvis.
- A flexible scoliosis can be managed by using the lateral supports in conjunction with the pelvic harness to provide 3-point positioning.
- The leg guides can accommodate up to 30° windsweeping, as well as abduction and adduction.
- Quickly and easily interfaces with a range of indoor and outdoor bases, including powered and non-powered wheelchairs, hi-low bases and strollers.



Available in 4 fabric colours

| Technical Details | | | |
|---------------------------------|--|--|--|
| | Size 1 | Size 2 | Size 3 |
| Age [approx] | 3-10 | 8-14 | 10-16 |
| User Weight | Min 18kg [40lbs] Max 50kg [110lbs] | Min 18kg [40lbs] Max 60kg [132lbs] | Min 18kg [40lbs] Max 70kg [154lbs] |
| User Height | Min 105cm [41"] Max 150cm [59"] | Min 127cm [50"] Max 168cm [66"] | Min 150cm [59"] Max 175cm [69"] |
| Seat Width | Min 200mm [8"] Max 325mm [13"] | Min 220mm [8.7"] Max 345mm [13.6"] | Min 250mm [9.8"] Max 400mm [15.7"] |
| Seat Depth | Min 270mm [10.6"] Max 420mm [16.5"] | Min 350mm [13.8"] Max 470mm [18.5"] | Min 350mm [13.8"] Max 470mm [18.5"] |
| Chest Width | Min 170mm [6.7"] Max 270mm [10.6"] | Min 170mm [6.7"] Max 270mm [10.6"] | Min 200mm [7.8"] Max 300mm [11.8"] |
| Backrest Height | Min 360mm [14.2"] Max 470mm [18.5"] | Min 460mm [18.1"] Max 570mm [22.4"] | Min 460mm [18.1"] Max 570mm [22.4"] |
| Backrest Angle | Prone 10° - Recline 25° | Prone 10° - Recline 25° | Prone 10° - Recline 25° |
| Lower Leg Length | Min 215mm [8.5"] Max 350mm [13.8"] | Min 315mm [12.4"] Max 470mm [18.5"] | Min 315mm [12.4"] Max 470mm [18.5"] |
| Armrest Height | Min 160mm [6.3"] Max 210mm [8.3"] | Min 210mm [8.3"] Max 260mm [10.2"] | Min 210mm [8.3"] Max 260mm [10.2"] |
| Tray Size | 550 x 480mm [21.6" x 18.9"] | 550 x 480mm [21.6" x 18.9"] | 550 x 480mm [21.6" x 18.9"] |
| Seat Height [footpedal chassis] | 340mm-655mm [13.4"-25.8"] | 340mm-655mm [13.4"-25.8"] | 340mm-655mm [13.4"-25.8"] |
| Seat Height [powered chassis] | 370mm-675mm [14.6"-26.6"] | 370mm-675mm [14.6"-26.6"] | 370mm-675mm [14.6"-26.6"] |

**Mygo Seat now
able to interface
with Quickie IRIS
Crash Tested Combination!**



LECKEY Everyday Activity Seat

Strong, durable, practical and adjustable seating system, providing moderate support for kids and young adults aged 1-18 years.



Product Features

- Ideal for everyday use- the soft-touch breathable fabric, well-padded ramped seat cushion, contoured back cushion and padded arm rests combine to offer the ultimate in comfort.
- The wide range of accessories and options make it easy to adjust the seat as the user grows and their needs change.
- The flip-up footplate, hi-low base [to floor], long armrests and 10° prone tilt-in-space support self-transfers.
- 4-point pelvic harness gives a secure, stable base of support.
- Height and angle adjustable chest laterals available.
- Easy care – covers can be removed and machine washed.



Available in 3 fabric colours and 2 vinyl versions

| Technical Details | | | | | | |
|--------------------------------|-------------------------|------------|-------------------------|------------|-------------------------|------------|
| | Size 1 | | Size 2 | | Size 3 | |
| | mm/kgs | inches/lbs | mm/kgs | inches/lbs | mm/kgs | inches/lbs |
| Max Load | 30kg | 66lbs | 50kg | 110lbs | 100kg | 220lbs |
| Weight - Seat & Hi-Low Chassis | 29.5kg | 65lbs | 34kg | 75lbs | 49kg | 108lbs |
| Seat Height [Hi-Low Chassis] | 250-600 | 9.8-23.6 | 320-600 | 12.6-23.6 | 450-600 | 17.7-23.6 |
| Seat Width | 160-280 | 6.3-11 | 260-380 | 10.2-15 | 340-500 | 13.4-19.7 |
| Seat Depth [Standard] | 150-280 | 5.9-11 | 230-380 | 9-15 | 330-480 | 13-18.9 |
| Seat Depth [Extended] | | | 300-450 | 11.8-17.7 | | |
| Lateral Distance - Chest | 150-220 | 5.9-8.7 | 170-240 | 6.7-9.4 | 240-370 | 9.4-14.6 |
| Back Height | 300-400 | 11.8-15.7 | 400-520 | 15.7-20.5 | 450-630 | 17.7-24.8 |
| Back Angle | 30° Supine 10° Prone | | 30° Supine 10° Prone | | 30° Supine 10° Prone | |
| Armrest Height | 150-250 | 5.9-9.8 | 200-300 | 7.9-11.8 | 180-325 | 7-12.7 |
| Lower Leg Length | 150-250 | 5.9-9.8 | 250-450 | 9.8-17.7 | 250-450 | 9.8-17.7 |
| Lower Leg Angle | -10° to +50° | | -10° to +50° | | -10° to +50° | |
| Footplate Angle | -15° to +10° | | -15° to +10° | | -15° to +10° | |
| Tilt-in-Space | 30° Supine 10° Prone | | 30° Supine 10° Prone | | 30° Supine 10° Prone | |

LECKEY Therapy Bench

Developed to support children in a neutral posture, this therapy bench with innovative activity cradles facilitates child active therapy by freeing up the hands.



Product Features

- Innovative, corseted activity cradle stabilises the pelvis in a neutral position, acting as an extra pair of hands for the therapist whilst giving the child the confidence to move, reach, and grasp.
- Height adjustable trunk support helps the child to sit upright with the head and spine stacked over the pelvis. This is the best and easiest way to learn sitting posture.
- Supportive, height adjustable therapist pad acts as extra support for the therapist who can rest against the rear support and gently encourage shoulder retraction.
- Facilitates lateral and anterior/posterior tilt to promote child activity.
- Height and depth adjustable bench making it versatile for both larger and smaller users.
- Acts as an aid to therapists in the Segmental Assessment of Trunk control [SATCo], providing a greater understanding of where to apply postural support to optimise function.



Therapy Bench

| Technical Details | |
|--------------------------------|--|
| One Size Only | |
| Max User Weight | 80kg (176.3lbs) |
| Age | 1-14 |
| Bench Footprint | 620mm x 740mm [63.7" x 29.1"] |
| Seat Base Height (from ground) | Min 330mm [12.9"] Max 460mm [18.1"] |
| Seat Base Width | 620mm [24.4"] |
| Seat Base Depth | 360mm [14.1"] |
| Footplate Height (from ground) | Min 20mm [0.7"] Max 110mm [4.3"] |
| Seat Base to Footplate Height | Min 220mm [8.6"] Max 440mm [17.3"] |

Activity Cradle

| Technical Details | | |
|-------------------|--|--------------------------------------|
| Cradle Size | Hip Circumference (mm) | ASIS - ASIS (mm) |
| Size 1 | Min 400mm [15.7"] Max 500mm [19.6"] | Under 120mm [4.7"] |
| Size 2 | Min 440mm [17.3"] Max 575mm [22.6"] | Min 130mm [5.1"] Max 160mm [6.2"] |
| Size 3 | Min 520mm [20.4"] Max 650mm [25.5"] | Min 160mm [6.2"] Max 180mm [7.0"] |
| Size 4 | Min 600mm [23.6"] Max 740mm [29.1"] | Min 190mm [7.4"] Max 200mm [7.8"] |
| Size 5 | Min 660mm [25.9"] Max 855mm [33.6"] | Over 200mm [7.8"] |

Standing and Walking

LECKEY Squiggles+

The Squiggles+ retains the best of the Squiggles stander but also offers heaps of added benefits so little ones can continue to enjoy standing and be part of the action.



Product Features

- Anatomically-correct abduction to stabilise the joint and provide optimum loading for early hip joint development. Squiggles+ offers 60° of hip abduction bilaterally in prone, upright or supine. It promotes proper alignment in the hips, knees and feet as pivot points are in line with the hips.
- Increased comfort and support with wrap-around knee supports. The comfortable medial knee supports distribute pressure and are especially suited to abducted standing. The padded wraps maintain the knee in neutral alignment.
- Providing greater confidence with transfers- the larger hip pads improve the carer's and child's confidence and feeling of security when transferring. The improved pivot chassis facilitates full horizontal loading.
- Activity tray with improved height and adjustability.
- Height-adjustable, foldable push handle.



| Technical Details | |
|--|--|
| | One Size Only |
| Age | 1-5 |
| Max User Weight | 22kg/48.4lbs |
| User Height | Min 730mm / 28.5 inches Max 1110mm / 43.5 inches |
| Chest Support | Angle adjustment +/- 30° Depth adjustment 50mm / 2 inches |
| Chest Support Height (footplate to armpit) | Min 560mm / 22 inches Max 940mm / 37 inches |
| Chest Width (distance between laterals) | Min 160mm / 6.3 inches Max 230mm / 9 inches |
| Hip Support (height from footplate to mid buttocks) | Min 300mm / 11.8 inches Max 530mm / 20.8 inches |
| Hip Width (distance between hip guides) | Min 160mm / 6.3 inches Max 230mm / 9 inches |
| Distance between Midline of Knees | Min 140mm / 5.5 inches Max 215mm / 8.4 inches |
| Footplate Angle | Plantarflexion 10° Dorsiflexion 10° |

Available in 3 fabric colours and a vinyl version



LECKEY Squiggles TT

Evolutionary early intervention standing frame that facilitates Targeted Training for children aged 1-5 years.



Product Features

- Maintaining the great features of the Squiggles+, the Squiggles TT goes a step further and facilitates Targeted Training. This evidence-based therapy is carried out with the child in an upright position, training one segment of their trunk at a time, allowing them to gain sequential control from the top down. Improved trunk control will help increase functional skills e.g. self-care and play. Better head control can help stop a child dribbling, improve their eye control as well as their ability to communicate.
- Squiggles TT is the only product specifically designed for Targeted Training therapy and which caters for all children. Horizontal loading is possible for complex children with poor head control and Squiggles TT is the only product to target all 7 trunk segments.
- Offers a fun and easy way for parents to carry out therapy at home. The snug supportive wraps are easy to use, the ratchet system with repeatable locking eliminates the fear of changing settings and the play ideas booklet provides support for home use.
- Dynamic base allows children to experience rocking movement whilst gaining upright postural control.
- Padded tray facilitates Targeted Training at the upper thoracic level.
- Conversion kit available for existing Squiggles+ users.

| Technical Details | |
|--|---|
| | One Size Only |
| Age [approx.] | 1-5 |
| Max User Weight | 22kg/48.4lbs |
| User Height | Min 730mm / 28.5 inches Max 1100mm / 43.3 inches |
| Trunk Support Height [Distance from footplate to armpit] | Min 160mm / 6.25 inches Max 820mm / 32 inches |
| Trunk / Pelvic Width | Min 135mm / 6.25 inches Max 250mm / 9 inches |
| Pelvic Support Height [Distance from footplate to mid-buttock] | Min 300mm / 11.8 inches Max 530mm / 20.8 inches |
| Distance between midline of knees | Min 140mm / 5.5 inches Max 215mm / 8.4 inches |
| Hip Width [distance between hip guides] | Min 160mm / 6.3 inches Max 230mm / 9 inches |
| Distance between Midline of Knees | Min 140mm / 5.5 inches Max 215mm / 8.4 inches |
| Footplate | Plantarflexion -10° Dorsiflexion -10° |
| Pivot chassis angle range | 90-180° |
| Tray Height | Min 480mm / 18.5 inches Max 825mm / 32.5 inches |



Available in 3 fabric colours and a vinyl version

LECKEY Squiggles Stander

Colourful, comfortable and supportive standing system for children aged 1-5 years.



Product Features

- 3-in-1 prone, supine and upright stander.
- Ideal for home use due to ease of adjustments.
- Minimalist child-focused design with fun bright pads blends easily into nursery/school environments.
- Low to the ground, encouraging socialisation and interaction at eye level.
- Individually adjustable chest and pelvic supports are cushioned for comfort and adjust easily in height, width and angle to suit individual needs, ensuring a comfortable and secure position.
- Frame can be easily disassembled for storage and transportation.



Available in 3 fabric colours and a Vinyl version

| Technical Details | |
|--|---|
| One Size Only | |
| Age | 1-5 |
| User Height | 730mm-1100mm [29.5"-43.3"] |
| Max. User Weight | 22kg (48lbs) |
| Chest Support Height (footplate to armpit) | 550mm-820mm [21.6"-32.2"] |
| Chest Width (distance between laterals) | 160mm-230mm [6.2"-9"] |
| Chest Support Adjustment | +/- 30° 50mm [1.9"] - Depth |
| Hip Support (height from footplate to midbuttocks) | 300mm-530mm [11.8"-20.8"] |
| Hip Width (distance between hip guides) | 160mm-230mm [6.2"-9"] |
| Distance between Midline of Knees | 140mm-215mm [5.5"-8.4"] |
| Footplate Angle | Plantarflexion 10° Dorsiflexion 10° |
| Tray Height (footplate to top of tray) | 555mm-825mm [21.8"-32.4"] - Prone 480mm-825mm [18.8"-32.4"] - Supine |
| Tray Angle Adjustment | 30° - Prone 40° - Supine |

LECKEY Mygo Stander

Strong, practical and adjustable standing system for kids aged 4-14 years.



Product Features

- 3-in-1 prone, supine and upright stander.
- Designed to improve upper limb function. The chest pad design replicates the shape of the rib cage which allows the shoulder girdle to freely protract and retract in prone, whilst in supine, the shoulder support promotes shoulder protraction. Furthermore, Mygo offers a lower tray height than most standers.
- Accommodates hip and knee flexion contractures of up to 25 degrees. The movement of the footplate and knee supports are linked which allows the user to stand with hips and knees flexed. This results in a more upright posture combined with the ability to support users whose contractures previously made them unsuitable for standing therapy.
- Use of the pelvic derotation belt puts the user's pelvis in a central and neutral position with minimum effort.



Available in 3 fabric colours

| Technical Details | | |
|---|--|--|
| | Size 1 | Size 2 |
| Age | 4-10 | 8-14 |
| User Height | 1080mm-1450mm [42.5"-57"] | 1350mm-1700mm [53.1"-66.9"] |
| Max. User Weight | 50kg [110lbs] | 60kg [132lbs] |
| Chest Support Height [footplate to armpit] | 760mm-1080mm [29.9"-42.5"] | 980mm-1300mm [38.5"-51.1"] |
| Chest Width [distance between laterals] | 160mm-280mm [6.2"-11"] | 200mm-320mm [7.8"-12.5"] |
| Hip Support Height [footplate to mid-buttocks] | 460mm-730mm [18.1"-28.7"] | 620mm-875mm [24.4"-34.4"] |
| Hip Width [distance between hip guides] | 220mm-330mm [8.6"-12.9"] | 280mm-390mm [11"-15.3"] |
| Footplate Angle | Plantarflexion 10° Dorsiflexion 10° | Plantarflexion 10° Dorsiflexion 10° |
| Max. Knee Contracture Accommodated | 25° | 25° |
| Tray Height Range [from top of hip pad] | 0mm-150mm [0"-5.9"] | 0mm-150mm [0"-5.9"] |
| Tray Angle Adjustment | 45°- Prone 45°- Supine | 45°- Prone 45°- Supine |

LECKEY Horizon

Robust, durable power stander for users aged 4 years to adult.



Product Features

- Accommodation of leg length discrepancy, scoliosis, lumbar lordosis, kyphosis, knee flexion contractures and asymmetrical head position. The frame's growth potential and the wide range of accessories accommodate changing need and potential asymmetry as the user grows.
- Perfectly sized for hoisting and supine standing. The large landing space for hoisted transfers coupled with support pads which provide a proximal fit to suit each user provides the perfect combination.
- Robust, stable and durable frame ensures ongoing confident product use.
- Individually adjustable chest, hip and knee pads as well as footplates provide unsurpassed levels of adjustability.
- Height adjustable hip laterals and pelvic belt, enabling de-rotation of the hips, alignment of the pelvis and centre of gravity over the base of support.



Available in 3 fabric colours

| Technical Details | | | |
|--|--------------------------------------|--------------------------------------|--------------------------------------|
| | Size 1 | Size 2 | Size 3 |
| Age | 4-10 | 8-14 | 14-adult |
| User Height | 1000mm-1400mm [39.3"-55.1"] | 1400mm-1800mm [55.1"-70.8"] | 1550mm-1800mm [61"-70.8"] |
| Max. User Weight | 50kg [110lbs] | 80kg [176lbs] | 100kg [220lbs] |
| Shoulder Height [footplate to top of pad] | 820mm-1150mm [32.2"-45.2"] | 1000mm-1450mm [39.3"-57"] | 1000mm-1450mm [39.3"-57"] |
| Hip Pad Height [footplate to top of pad] | 450mm-750mm [17.7"-29.5"] | 750mm- 950mm [29.5"-37.4"] | 750mm-1150mm [29.5"-45.2"] |
| Hip Pad Width | 175mm-270mm [6.8"-10.6"] | 210mm-390mm [8.2"-15.3"] | 310mm-470mm [12.2"-18.5"] |
| Knee Height (footplate to middle of pad) | 250mm- 450mm [9.8"-17.7"] | 400mm-550mm [15.7"-21.6"] | 400mm-650mm [15.7"-25.5"] |
| Width Between Chest Laterals | 200mm-300mm [7.8"-11.8"] | 240mm-380mm [9.4"-14.9"] | 340mm-460mm [13.3"-18.1"] |
| Footplate Angle Adjustment +/- | 10° | 10° | 10° |
| Angle Adjustment | 0°-90°- Prone 0°-90°- Supine | 0°-90°- Prone 0°-90°- Supine | 0°-90°- Prone 0°-90°- Supine |
| Tray Angle Adjustment | 0°-40°- Prone +/- 20°/30°- Supine | 0°-40°- Prone +/- 20°/30°- Supine | 0°-40°- Prone +/- 20°/30°- Supine |

LECKEY Pronestander

Robust, adjustable prone standing system for users aged 1-18 years.

Product Features

- Flip-up sandals support users completing a standing transfer directly into the product.
- Long base ensures maximum stability for even the most active users, providing caregivers with peace of mind.
- Unique de-rotation belt and pelvic pad work together to keep the user in a central position within the frame with equal pressure over each foot.
- Individually adjustable chest, hip, knee and sandal supports to meet a wide range of clinical needs.
- Easy to set up and adjust with five incremental secure prone positions from 85 to 45 degrees.



Size 1 only available in Meadow Fabric
Size 2 & 3 only available in Mondeo Fabric

| Technical Details | | | |
|------------------------------|----------------------------|---------------------------|-----------------------------|
| | Size 1 | Size 2 | Size 3 |
| Age | 1-5 | 4-10 | 9-18 |
| User Standing Height | 800mm-1151mm [31.4"-45.3"] | 1100mm-1448mm [43.3"-57"] | 1400mm-1849mm [55.1"-72.7"] |
| Max. User Weight | 60kg [132lbs] | 80kg [176lbs] | 90kg [198lbs] |
| Chest Pad Height | 470mm-680mm [18.5"-26.7"] | 680mm-975mm [26.7"-38.3"] | 975mm-1350mm [38.3"-53.1"] |
| Prone Angles | 85°, 75°, 65°, 55°, 45° | | |
| Width between Chest Laterals | 152mm-254mm [5.9"-10"] | 203mm-305mm [7.9"-12"] | 254mm-356mm [10"-14"] |
| Tray Sizes | 540mmx400mm [21.2"x15.7"] | 590mmx425mm [23.2"x16.7"] | 644mmx470mm [25.3"x18.5"] |

LECKEY Freestander

Robust, upright standing system for users aged 1-18 years.

Product Features

- Flexibility to alter the frame to meet individual user needs e.g. pelvic support can be removed for users with trunk control, conversely, the tray height and angle can be increased where more support is required. A simple, active set-up is also possible for more able users e.g. hemiplegic or GMFCS II/III.
- Straightforward, intuitive adjustment as all supports run off two vertical uprights.
- Pelvic belt achieves a unique position which enables de-rotation of the hips, alignment of the pelvis and centre of gravity over the base of support as well as accommodation or correction of flexion at the hips and knees.



Size 1 only available in Meadow Fabric. Size 2 & 3 only available in Mondeo Fabric.

| Technical Details | | | |
|------------------------------|----------------------------|---------------------------|-----------------------------|
| | Size 1 | Size 2 | Size 3 |
| Age | 1-5 | 4-10 | 9-18 |
| User Standing Height | 800mm-1151mm [31.4"-45.3"] | 1100mm-1448mm [43.3"-57"] | 1400mm-1849mm [55.1"-72.7"] |
| Max. User Weight | 60kg [132lbs] | 80kg [176lbs] | 90kg [198lbs] |
| Chest Pad Height | 470mm-680mm [18.5"-26.7"] | 680mm-975mm [26.7"-38.3"] | 975mm-1350mm [38.3"-53.1"] |
| Width between Chest Laterals | 152mm-254mm [5.9"-10"] | 203mm-305mm [7.9"-12"] | 254mm-356mm [10"-14"] |
| Tray Sizes | 540mmx400mm [21.2"x15.7"] | 590mmx425mm [23.2"x16.7"] | 644mmx470mm [25.3"x18.5"] |

LECKEY Todstander

Simple and easy to use, upright standing system for kids aged 1-5 years.

Product Features

- Can be used as a stander or by therapists for standing assessments. Promotes active standing, movement and weight shift.
- Compact and easily transported to different environments e.g. from therapy centres to the home.
- Provides a gentle introduction to rehab equipment via a simple, easy-to-use, child-friendly design.

| Technical Details | |
|-------------------|----------------------------|
| One Size Only | |
| Age | 1-5 |
| Chest Band Height | 470mm- 640mm [18.5"-25.1"] |
| Base Depth | 450mm [17.7"] |
| Base Width | 450mm [17.7"] |
| Max. User Weight | 20kg [44lbs] |



Available in 1 fabric colour combination

LECKEY MyWay

MyWay is the leader of the pack when it comes to maximising opportunities for independent movement and giving children the freedom to explore their world.



Product Features

- The unique corset-like harness gives core trunk stability for optimal stepping pattern and active movement. The MyWay harness provides the support required for users to develop head and limb control. Its hug-like effect increases sensory feedback, dampens involuntary movement and increases confidence and the chance of succeeding.
- MyWay's open-front, hands-free frame maximises opportunities for exploration and participation. The open-design also allows therapists to get close to the child and work with their lower limbs whilst they are in MyWay, facilitating guided movement.
- The hoistable harness makes transfers easier- direct transfers eliminate the need for additional hoisting slings and make the most out of therapy time. Hoisting in the MyWay harness also provides superior comfort to the user as it does not slip into the armpits.
- Multi-functional castors, lightweight frame and prone, gross and fine height adjustment combine to help initiate stepping and control movement.
- Range of optional accessories available to meet individual requirements, including headrest, ankle guides and saddle.



Each harness size is a different color

| Technical Details | | | |
|---------------------------------|--|--------------|--------------|
| | Size 1 Frame | Size 2 Frame | Size 3 Frame |
| Age (approx.) | 1-5 | 4-11 | 10-16 |
| User Shoulder Height (to floor) | 635-920mm | 870-1190mm | 1070-1490mm |
| User Saddle Height (inside leg) | 220-440mm | 400-610mm | 550-810mm |
| Max. User Weight | 25kg | 50kg | 80kg |
| Frame Weight | 12 [Size 1 frame and upper support unit do not separate] | 6.5kg | 9.5kg |
| Upper Support Unit Weight | | 10.5kg | 11kg |
| Frame Width (mm) | 570mm | 670mm | 720mm |
| Frame Length (mm) | 720mm | 850mm | 1100mm |

| Harness | | | | | | |
|---------------------------|--------|--------|--------|--------|--------|----------|
| | Size 1 | Size 2 | Size 3 | Size 4 | Size 5 | Size 6 |
| Age (approx.) | 1-4 | 2-6 | 5-9 | 8-12 | 12-15 | Up to 16 |
| Waist at Navel Level (cm) | 46-56 | 50-60 | 55-65 | 60-70 | 70-80 | 80-90 |
| Colour | Purple | Green | Blue | Red | Orange | Grey |

LECKEY MyWay Pedal

MyWay Pedal is a revolutionary elliptical trainer which helps improve strength and functional ability. Attaching to the size 2 and 3 MyWay frame, MyWay Pedal is ideal for giving on-the-go kids the opportunity to be active, have fun and join in with their friends.

Product Features

- Targets and strengthens anti-gravity muscles to help reduce the risk of hip and knee contractures and improve function. Pedalling strengthens anti-gravity muscles (calf, quads and hamstring)- into extension unlike a trike which encourages a flexed position. MyWay Pedal ensures the quads are worked exactly where strength is often lacking, at the end of range. Working the end of range will reduce the risk of hip and knee contractures which may limit standing and walking ability.
- Facilitates elliptical action more like stepping to help support gait therapy goals. MyWay Pedal provides a sensori-motor experience similar to gait. Unlike the circular motion of cycling, the elliptical pedalling of MyWay Pedal gives an extended stretch of the lower limbs as it involves a longer forwards/backwards movement rather than an up/down movement. It helps improve loading around the hip to aid joint development.
- Use in conjunction with the Leckey Connect App to chart performance and motivate the child.



The app provides a valuable way of measuring performance and progress on a child's speed, distance and duration of elliptical pedalling. It is available on Android Smartphone, tablet and the Leckey online portal.

- Attaches to the size 2 and 3 MyWay frame and can be fitted with the user in or out of the frame.
- Two sandal sizes available.
- Height adjustable and removable handlebar.



Other Products

LECKEY Advanced Bath Chair

The Leckey Advanced Bath Chair provides postural support while ensuring the individual is safe and secure during bathing.

Product Features

- Provides the user with appropriate levels of postural support ensuring they are safe and secure at all times during bathing.
- Lightweight and robust frame.
- Simple and easy-to-use adjustments mean that the product can be set up in a variety of positions in seconds.
- Height adjustable legs.
- Breathable cover, removable and machine washable.
- Back angle adjustment from 0-90 degrees.
- Adjustable hip belt.
- Chest belt with integrated lateral supports.
- Easy to use hand adjustments.
- Folds flat for storage.



| Technical Details | | | | |
|-------------------|---------------|---------------|---------------|---------------|
| | Size 1 | Size 2 | Size 3 | Size 4 |
| Age | 1-5 | 4-9 | 8-14 | 12-18 |
| Min User Height | 750mm [30"] | 950mm [37"] | 1150mm [45"] | 1350mm [53"] |
| Max User Height | 1050mm [41"] | 1250mm [49"] | 1450mm [57"] | 1650mm [65"] |
| Max User Weight | 72kg [159lbs] | 72kg [159lbs] | 72kg [159lbs] | 72kg [159lbs] |



Available in 2 fabric colours

LECKEY Shower Trolley

The Shower Trolley complements the Advanced Bath Chair which can be converted within seconds into a shower chair.

Product Features

- Converts the Advance Bath Chair to a versatile shower chair in seconds.
- Manufactured from medical grade stainless steel.
- Robust frame and high quality lockable castors.
- Can be easily and safely manoeuvred over wet and slippery surfaces.
- Can be adjusted to fit all four sizes of the Leckey Advance Bath Chair.



| Technical Details | |
|------------------------|------------------------|
| | One Size Only |
| Age | 1-18 |
| Length | 900mm-1150mm [36"-45"] |
| Width | 600mm [24"] |
| Seat Height from Floor | 650mm-700mm [26"-28"] |

LECKEY Sleepform

Strong, durable, practical and adjustable sleeping system for kids aged 0-18 years.

Product Features

- Close proximal fit cocoons the child and makes them feel safe and secure. Adapts quickly and easily to changing body shape.
- Compact and portable for visits to family, going on holiday etc. Suitable for use on a variety of surfaces.
- Supports 24-hour postural management as can be used during the day to support changes of position.
- Can be used in cots and single beds and suits all sleeping positions.
- Air flow mattress.
- Grid sheet to record best sleeping positioning.
- Temperature control sheet available.

| Technical Details | | | | |
|-------------------|--------|--------|------------|------------|
| | Size 1 | Size 2 | Size 3 | Size 4 |
| Age | 0-1 | 1-5 | 1-5 | 5-18 |
| Suitable for | Cot | Cot | Single Bed | Single Bed |



Supporting Children





©2021 Sunrise Medical Leckey Range AU Rev D

Sunrise Medical Australia Pty Ltd
11 Daniel Street, Wetherill Park, NSW 2164
Phone: (02) 9678 6600
Email: enquiries@sunrisemedical.com.au
Website: www.SunriseMedical.com.au