

Intelligent Rotation In Space Technology

QUICKIE[®]
IRIS





Care Giver Options

The IRIS tilt wheelchair is available with numerous options to make your job as caregiver that much easier. For cable-free tilting, choose the Foot Release Tilt Actuator. The hands-free Attendant Wheel Lock also uses a foot pedal to engage. Standard with every IRIS, the Angle Adjustable Stroller Handle can be easily adjusted to differing caregiver heights for easy handling.



Streamlined, Lightweight Design

- The IRIS is one of the lightest positioning wheelchairs on the market and is even easier to transport and maneuver.
- Designed with a low caster fork and small caster wheel options, the IRIS can achieve a seat-to-floor height as low as 12.5".
- The IRIS's incredibly low seat-to-floor height provides easier access to tables, desks and counters as well as enhanced floor access for foot propulsion.



MONO™ Backrest System

The MONO Backrest System's single-post offers limitless seating and positioning capabilities. By mounting the backrest to a single central post, the backrest and chair widths may be selected independently and hardware interference is eliminated. The MONO Backrest System offers an optional dynamic function, 35° angle adjustment, and folds forward for easy transport.



iXpress Power Add On Compatible

The iXpress offers a high performing, reliable Power Add On System that provides the convenience & transportability of a manual wheelchair but with speed of a power chair.

The ultimate tilt-in-space wheelchair, the **Quickie IRIS** features up to 55° of Intelligent Rotation in Space Technology, the broadest range of positioning and caregiver options, numerous adjustments, and it's one of the lightest tilt-in-space wheelchairs available today.

Intelligent Rotation in Space (IRIS)

Intelligent Rotation in Space technology means the seat frame rotates around the user's center of gravity. This results in a smooth tilting motion that is easy for you to control and allows for the shortest possible wheelbase that is easy for you to maneuver.

The IRIS manual wheelchair features 40° or 55° rotation ranges that can be set to achieve up to 45° or 60° of posterior rotation for superior positioning.

True Fit Program

If the need arises to significantly modify the chair, the Quickie Iris comes with a True Fit program which provides a one time free part offer within 5 years of original chair ownership.

Fold Down Back
Canes to increase
portability

Up to 55° of rotation



Intuitive Hanger design -
easy to use and remove.
Swing away and Elevating
will interchange

Rocker Arm designed
for smooth transitions



Easy to use Foot
Release Tilt Actuator

Wide Range of Frame
Colours Available



Z-finity footrest with
infinite adjustment in
both height and angle



Jay Lower Extremity
Support Option

QUICKIE[®] IRIS

QUICKIE[®] IRIS

SPECIFICATIONS

Frame Style

Tilt-in-Space

Transport Weight

13 kgs

Overall Weight

17.5 kgs

User Weight Capacity*

136 kgs Standard

182 kgs Heavy Duty

Seat Widths

14" – 22"

Frame Depths

14" – 22"

Front Seat-To-Floor Height

12.5" – 19.5"

Tilt Range

-5° to 35°, 0° to 40°, and +5° to 45°

OR

-5° to 50°, 0° to 55°, and +5° to 60°



Sunrise Medical Australia Pty Ltd
6 Healey Circuit, Huntingwood, NSW 2148
Phone: (02) 9678 6600
Fax: (02) 9678 6655
Email: enquiries@sunrisemedical.com.au
www.sunrisemedical.com.au

© Sunrise Medical. Quickie Iris Rev D.AU