

A black and silver Permobil 2GX wheelchair is shown from a side profile, parked on a paved surface. The wheelchair features a folding frame, large rear wheels with green spokes, and smaller front wheels. The background is a blurred outdoor setting with a building and greenery. The image is overlaid with a large, abstract geometric graphic in shades of teal and blue on the right side.

permobil

# 2GX

## Titanium performance.

The 2GX combines the advantages of a modular folding frame with the high strength and vibration dampening characteristics of titanium. Add full adjustability and the 2GX is designed to empower - no matter where you choose to roll.



## Technical Specifications

Frame Features	Titanium Box Frame Folding Modular Design Full Adjustability
Total Width Folded*	10.75"
Total Chair Width*	23"
Seat Width	12" – 22"
Seat Depth	12" – 20"
Front Seat Height	13.5" – 21"
Rear Seat Height	13" – 20"
Front Frame Angle	Fixed Front - 70°, 80°, or 85° Swing Away - 70°, 80°, 85°, or 90°
Seat Back Type	<i>Folding Seat Back</i> Adj. Height 11" – 20.5" Adj. Angle 80° – 101°  <i>Fixed Seat Back</i> Adj. Height 11" – 20"

Camber	0°, 2°, 4°, or 6°
Centre of Gravity Adjustment	6"
Transit Options	Transport Bracket
Transport Weight*	Fixed Front - 5.9 kg Swing Away - 7 kg
Weight Limit	120.2 kg Standard 158.7 kg Heavy Duty
Warranty	Five Years on Frame One Year on TiLite Components

\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.

## Modular Design

Optimise your ride with custom configuration. This chair's modular design allows you to add options including after market seating systems, switching fixed and swing away front ends, and growing chair width and / or depth to meet changing user needs.

## Custom Configuration

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.


## Colour Options

### High Gloss Paint



TiLite Red    Burgundy    Neon Pink    Super Wet Black

### Matte Paint



Deep Blue Ice    Titanium Grey    Black Lava

### Metallic Paint



Salsa Red Metallic    Tangerine Metallic    Acid Green Metallic    Ocean Blue Metallic

### Candy Paint



Ruby Red    Azure Blue    Electric Plum

### Natural



Satin Titanium

### Pearlescent Paint



Midnight Blue Pearl    White Pearl