



## Introducing the Casero bed series

*Meet the bed designed for your needs and your home*



*Delivering home style luxury beds  
across all care settings*







## Features of the future

The Casero bed series showcases sophisticated, stylish and intelligent solutions to suit both client's needs and their homes.

Developed especially for homecare and aged care environments, the Casero series offers 4 different bed styles covering aged care and community clients, alongside additional and highly functional features to reduce falls risks as well as options for bariatric users and couples.

The Casero series is modular and configurable to your client's requirements and includes classic and contemporary designs, upholstery choices and high quality features.

The European designed 4-part telescopic side rail option offers versatility, support, safety, protection – and overall peace of mind. With an extensive range of options and accessories, the Casero series truly reimagines and merges perfectly the functional needs of our clients and the style that is important to them.







# Casero styles & finishes

The Casero range is modular and configurable to your client's requirements and includes classic and contemporary designs, upholstery choices and high quality features.

Headboard styles		Colour options	
The Oxford		Oatmeal	
		Granite	
The Belgrave		Slate	
		Ebony (vinyl)	

\*All coverings available in each headboard style

# Casero features & accessories

## The Duo

Width and length adjustable (Single or King Single sizes)  
Height range 184mm - 765mm  
Max user weight 250kg  
(BD6640)

---

## The Mode

Available in Single and King Single sizes  
Disassembles for easy transportation ideal for homes with restricted accessibility  
Low height range ideal for falls risk clients  
Height range 110mm - 740mm  
Max user weight 200kg  
Single (BD5670)  
King Single (BD5680)

---

## The Vogue

Available in Single and King Single sizes  
Low height range ideal for falls risk clients  
Height range 110mm - 740mm  
Max user weight 250kg  
Single (BD5210)  
King Single (BD5220)

---

## The Plus

Width and length adjustable (Double or Queen sizes)  
Suitable for couples or clients requiring greater sleeping area or increased weight capacity  
Height range 250mm - 660mm  
Max user weight 450kg  
(BD5290)

## Features

Electric Hi/Low Functions  
Electric Back Elevation  
Electric Knee Break  
Electric Trendelenburg\*  
Built-in 200mm length extension

## Accessories

Self-help Pole (BZ0325)  
Telescopic 4-Part Side Rails (BZ0443)  
Length Horizontal Side Rails (BZ0302)  
Drop-down Quarter Rails (BA5680)  
Battery Back-up System (BZ0335)\*  
Nurse Control Panel (BZ0337)\*  
Under-Bed Light (BZ0336)\*

## Warranties

Steel frame - 10 years  
Electrics - 5 years  
Hand control - 12 months  
Upholstered surround incl head and foot boards - 12 months  
\* not available on The Plus model



# Quietly caring

The Casero bed series showcases sophisticated, stylish and intelligent solutions to suit both client's needs and their homes

---

## Australia wide locations

Sydney | Newcastle | Brisbane  
Melbourne | Canberra | Griffith

02 9714 1111  
info@carequip.com.au  
**carequip.com.au**

Head office  
5 Millennium Court  
Silverwater NSW 2128

---

