


# [todo] DRIVE

START A NEW ADVENTURE!



*Todo Works creates appropriate technology to resolve the issues faced by people with reduced mobility. Each person's need is just as important as everyone's need. If there is something only one person needs, that is absolutely what we have to do. We continue working for what is needed for an individual as well as for all.*

[todo]  
WORKS





- 

### ***WORLD'S EASIEST***

Have freedom of mobility by mounting Todo Drive which is composed of Motor, Controller, Main Unit, and Battery. **A 5mm hex key is only necessary tool** for you to install it. Furthermore, it simply transforms manual and power mode by pushing and pull levers.



FIRST  
ADVENTURE

**SO  
EASY**





••

## ***WORLD'S LIGHTEST***

Feel comfortable and convenient like a power wheelchair **by adding this 5kg kit on your wheelchair**. It can be easily loaded in small compact cars because the folding function of the existing wheelchairs is unrestricted.



SECOND  
ADVENTURE

**JUST  
5kg**





...

## ***WORLD'S MOST UNIVERSAL***

There are about 1,500 kinds of wheelchairs in the world. The universal design of **Todo Drive is fitted to most of these wheelchairs.** Todo Drive is simply detachable and attachable. No problem even if you change your wheelchair.



THIRD  
ADVENTURE


**FITS  
ALL**



....

## ***SMART IoT READY***

There have various IoT sensors in Todo Drive. These sensors make it easy to use additional functions such as **remote control and personal setting via your phone**. The new functions of Todo App will be regularly updated in Appstore.





FOURTH  
ADVENTURE

# SMART IoT





***SIMPLY FIVE***

Only 5parts needed to get you stated.



FIFTH  
ADVENTURE

# SIMPLY FIVE



.....

## ***SPECIFICATION***

**RANGE** 7km( $\pm 5\%$ )\*

**CLIMBING CAPABILITY** Up to 6\*\*

**SPEED** 0.2~6km/h\*

**RECOMMENDED USER WEIGHT** Up to 90kg

**OPERATION TEMPERATURE** -10°C~60°C

**MOTOR RATING** 2 X 180W

**TYPE OF BATTERY** Lithium-ion Rechargeable Battery

**BATTERY PACK** 14.4V, 10Ah, 144Wh

**BATTERY WEIGHT** 0.99kg

**MAIN UNIT WEIGHT** 0.42kg

**MOTOR WEIGHT** 2 x 1.4kg

**TOTAL WEIGHT** 5kg\*\*

**BLUETOOTH SPEC** 2.4GHz Bluetooth 4.0

\*May differ depending on weight, environment, usage pattern, etc.

\*\*May change depending on the option.





SIXTH  
ADVENTURE



.....

## **START A NEW ADVENTURE!**

*Todo Drive is a power add-on for manual wheelchair to expand the range of mobility for wheelchair users by adding a simple kit to frame.*

**Start a new adventure, new experience,  
new every day with Todo Drive.**