



Hoist Accessories

Jordan Frame (& Strap Kit)

- Supine lifting device including 2 piece frame, gliders & strap kit for hoist
- 2 sizes – wide (SWL 200kg) and extra wide (SWL 230kg)
- Suitable for yoke hook attachment (preferably 4 Point Yoke)

Available on: Alto, Kratos and other large non-Allegro hoists

Indications: Lifting patients requiring spinal or orthopaedic immobility

Useful for: Floor lifts where spinal or orthopaedic injury is suspected

Weigh Scale

- Universal weighing module with tare function and kg/lb conversion
- SWL up to 300kg
- Suitable for yoke or pivot attachments

Available on: All Allegro & Allegro Lite hoists plus some non-Allegro hoists

Indications: Where weighing of patients is difficult due to the level of support required

Useful for: Patients requiring constant monitoring of weight for general wellbeing

Jordan Frame



Weigh Scale



**Innovation in
Patient Handling**

Allegro
CONCEPTS PTY LTD

Hoist Accessories

Allegro Concepts offer a broad range of accessories that interface with the Allegro steel frame lifting hoists, the Allegro Lite aluminium frame lifting hoists plus many non-Allegro lifting hoists.

Some accessories [such as the Jordan Frame] are designed to increase the functionality of existing patient lifting equipment. In this case a regular patient lifting hoist with yoke hook attachment can be used to lift patients requiring total spinal immobility e.g. unstable quadriplegia immediately post trauma.

Other accessories [such as the Sling Bag] help to make the process of lifting and transferring patients a more ergonomic task. Time is saved with the carer having ready access to the lifting sling at all times.

Allegro Concepts Hoist Accessories are designed and manufactured in Australia, with a view to providing healthcare facilities with customized patient transfer packages that meet the needs of each individual client.

Sling Bag

- Heavy duty vinyl bag for storing slings and other patient handling accessories (e.g. Swift Slide Sheets)
- Attaches to the push handles of the hoist with 2 Velcro tabs
- Custom sizes and monogramming available on request

Available on: All Allegro and non-Allegro patient lifting hoists

Indications: Where storage of and immediate access to the sling is required

Useful for: Storing slings, Swift Slide Sheets and Patient Transfer Belts on the end of patients beds

Tempo Handle Kit

- Pair of bolt on handles with non slip rubber grip
- Designed to provide optimum grip position on Salsa standing transfer hoist
- Multi adjustable in angle, height and distance from patient

Available on: Allegro Salsa standing transfer hoist

Indications: Where the regular grip position is not comfortable or ergonomically placed for the patient

Useful for: Hemiplegic patients who require optimum gripping position with the one functional hand

(Salsa) Leg Strap Knob Kit

- Hardware for the attachment of leg straps when custom fitted to an Allegro standing sling
- Fitted to upper clevis of lifting actuator (can be retrofitted)
- Allows the use of the Allegro Salsa for non-standing transfers

Available on: Allegro Salsa standing transfer hoist

Indications: Patients with sitting balance but limited (or fluctuating) weight-bearing capacity

Useful for: Conditions such as Multiple Sclerosis where transfer ability can fluctuate within a short space of time

Leg Protector

- Clear Polycarbonate plastic covers for protecting paint finish
- Fitted and removed without tools
- Easily cleaned for infection control

Available on: Allegro steel frame patient lifting hoists

Indications: Facilities which rely heavily on patient lifting hoists and are concerned with general upkeep of equipment

Useful for: Preventing rusting of patient lifting hoists used in wet environments

Sling Bag



Tempo Handle Kit



Knob Kit



Leg Protector



Allegro Concepts Pty Ltd

30-32 James Street Lidcombe NSW 2141

Phone: (02) 9749 7812 Fax (02) 9749 2144

Email: sales@allegroconcepts.com.au

Allegro
CONCEPTS PTY LTD
www.allegroconcepts.com.au