

**LECKEY**



# Squiggles+

Perfect blend of style and  
function for little ones

Supporting Children

# Building on the best of the much-loved Squiggles stander

The incredibly versatile **three-in-one Squiggles stander**, with its minimalist child-focused design, has always been a standout when it comes to happy parents and happy children.

Noted for being **low to the ground**, the Squiggles stander is ideal for encouraging socialisation and interaction at eye level, while offering **unsurpassed support, adjustability and easy storage**.

Retaining the best of the Squiggles stander, we have now created the **Squiggles+** with heaps of added benefits so little ones can continue to enjoy standing and be part of the action.



## Squiggles+

Perfect blend of style and function for little ones

## 1 ANATOMICALLY-CORRECT ABDUCTION FOR EARLY HIP JOINT DEVELOPMENT

- Stabilises the joint and provides optimum loading
- Offering 60° of hip abduction bilaterally in prone, upright or supine
- Promoting proper alignment in the hips, knees and feet as pivot points are in line with the hips
- Hip laterals that shift in sync with the child's leg allowing the hips to abduct and preventing discomfort

## 2 WRAP-AROUND KNEE SUPPORTS FOR EXTRA COMFORT

- Medial knee supports to distribute pressure, especially suited to abducted standing
- Padded wraps to maintain knee in neutral alignment
- Can be rotated to provide extra lateral support to the knee for low tone children

## 3 TRANSFER WITH CONFIDENCE

- Larger hip pads to improve carer's and child's confidence and feeling of security when transferring
- Improved pivot chassis to facilitate full horizontal loading





# Technical Information

## Stander Dimensions

Age (approx.)	1-5
Max User Weight	22kg/48.4lbs
User Height	Min 750mm/29.5 inches   Max 1110mm/43.5 inches
Chest Support	Angle adjustment +/- 30° Depth adjustment 50mm/2 inches
Chest Support Height (footplate to armpit)	Min 550mm/21.5 inches   Max 820mm/32 inches
Chest Width (distance between laterals)	Min 160mm/6.3 inches   Max 230mm/9 inches
Hip Support (height from footplate to mid buttocks)	Min 300mm/11.8 inches   Max 530mm/20.8 inches
Distance between midline of Knees	Min 140mm/5.5 inches   Max 215mm/8.4 inches
Footplate Angle	Planterflexion 10° / Dorsiflexion 10°

## Tray Size

Tray Height - Prone (footplate to top of tray)	Min 555mm/21.9 inches   Max 800mm/31.5 inches
Tray Height - Supine (footplate to top of tray)	Min 480mm/18.5 inches   Max 800mm/31.5 inches
Tray Angle Adjustment	Prone 30° / Supine 40°

## Support Frame & Pivot Chassis Dimensions

Stander support frame	Weight 7.5kg / 16.5lbs   Width 390mm / 15.4 inches Length 680mm / 26.8 inches   Height 390mm / 15.4 inches
Chassis footprint (unfolded)	Weight 6.5kg / 14.3lbs   Width 590mm / 23.2 inches Length 950mm / 37.4 inches   Height 515mm / 20.3 inches
Chassis footprint (folded)	Weight 6.5kg / 14.3lbs   Width 590mm / 23.2 inches Length 660mm / 26.0 inches   Height 515mm / 20.3 inches
Chassis angle range	90 - 180°   90-70°

Leckey  
19c Ballinderry Road  
Lisburn BT29 2SA  
Northern Ireland

[+44] 28 9260 0750  
hello@leckey.com

leckey.com



**LECKEY**