

# NIMBO

## Lightweight Gait Trainer



The Nimbo Posture Walker supports the user from the rear encouraging upright posture and trunk extension, suitable for children with Cerebral Palsy and other neurological conditions.

- Supports user from rear encouraging upright posture
- Sturdy, yet lightweight aluminium frame
- Tool free height adjustment to suit various users
- Latex free handles with flanged ends

[www.drive-medical.com.au](http://www.drive-medical.com.au)

**drive**

**DeVilbiss**  
HEALTHCARE

# NIMBO

## Lightweight Gait Trainer



### Product features

- Supports user from rear encouraging upright posture and trunk extension - suitable for children with Cerebral Palsy and other neurological conditions
- Sturdy, yet lightweight aluminium frame is available in attractive bright colours
- Tool free height adjustment to suit various users
- Latex free handles with flanged ends allow comfortable, firm handgrip
- Soft rubber wheels adhere to any surface and prevent Nimbo from sliding backwards
- One directional rear wheels with aluminium ratchet and pin mechanism.
- Front swivel wheels enhance manoeuvrability and ease of turning
- Anti-reverse override bracket disengages the one directional rear wheels allowing forward & reverse mobility
- Locking mechanism easily switches front wheels from swivel to non-swivel for greater control
- Each wheel is independently adjusted for flexibility
- Skid-resistant, non-marking and easy rolling wheels ease manoeuvrability and increase stability
- Folds easily for convenient storage and transportation
- Available in five sizes with a range of optional accessories
- 2 year limited warranty

### Accessories



#### 1. FOREARM PLATFORM

Constructed of lightweight aluminium and height/depth adjustable. The handgrips extend forward or back to the position needed, supinate, pronate and are angle adjustable. Each side is individually adjustable to allow custom positioning. The platforms Velcro straps prevent the arms from slipping. They can be mounted either towards the inside or outside of the walker for lateral adjustment. Spacers are included for additional lateral range.

CODES: KA1035FPS (SMALL)  
KA1035FPL (LARGE)



#### 2. PELVIC STABILIZER

The pelvic stabilizer is adjustable in width/height and features back & lateral pads and a support belt. It is recommended for users who have an asymmetrical posture when standing or walking and therefore, cannot stay centered in the walker.

CODE: KA1050



#### 3. SEAT & SEAT HARNESS

Seat mounts on standard walker. They fold up for standing and walking and fold down for convenience.

Adjustable sling seat supports the user that cannot stand completely on their own.

CODES: KA1285 (SEAT FOR KA1200N)  
KA2285 (SEAT FOR KA2200N)  
KA3285 (SEAT FOR KA3200N)  
KA4285 (SEAT FOR KA4200N)  
KA5285 (SEAT FOR KA5200N)  
CE1070S (SEAT HARNESS FOR KA1200, KA2200 AND KA3200)  
CE1070L (SEAT HARNESS FOR KA4200 AND KA5200N)



#### 4. 9" OUTDOOR WHEELS

CODES: KA8100 (PAIR FOR KA2200 - KA5200)  
KA8100S (PAIR FOR KA1200)

#### 5. DIAL-A-SPEED TABS

Provides variable resistance on one or both rear wheels. They allow user to control gait trainer speed and aligns asymmetrical gait patterns.

CODE: KA1400 (PAIR)



#### 6. BUMPER GUARDS

CODE: KA8000 (PAIR)

CODE	DESCRIPTION	HANDLEBAR HEIGHT RANGE	WIDTH OPENED AT BASE	WIDTH INSIDE HANDLE GRIP	DEPTH OPENED AT BASE	WEIGHT	WEIGHT CAPACITY
KA1200N	Posture Walker - Tyke, Goldenrod Yellow	39 - 50cm (15.5" - 19.5")	61cm (24")	34cm (13.5")	51cm (20")	4kg (9lb)	34kg (5st)
KA2200N	Posture Walker - Junior, Lime Green	48 - 64cm (19" - 25")	61cm (24")	34cm (13.5")	69cm (27")	5kg (10lb)	39kg (6st)
KA3200N	Posture Walker - Youth, Cornflower Blue	50 - 77cm (23" - 30.5")	66cm (26")	38cm (15")	79cm (31")	6kg (13lb)	68kg (11st)
KA4200N	Posture Walker - Young Adult, Midnight Blue	71 - 91cm (28" - 36")	67cm (26.5")	42cm (16.5")	86cm (34")	6kg (14lb)	86kg (14st)
KA5200N	Posture Walker - Adult, Flame Red	92 - 104cm (36" - 41")	74cm (29")	46cm (18")	94cm (37")	7kg (16lb)	113kg (18st)