

Room to grow

Small chair with mini kit at its smallest dimensions

- 1 Backrest height 12½"
- 2 Trunk support width 5½"
- 3 Seat depth 7"
- 4 Seat width 7"
- 5 Armrest height above seat 5"
- 6 Seat height above footboard 7"




Small chair at its largest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 11½"
- 3 Seat depth 12"
- 4 Seat width 12"
- 5 Armrest height above seat 7½"
- 6 Seat height above footboard 12"




Small standard base
R820 Rifton Activity Chair

User dimensions (inches)	
Height with mini kit	32–48 28–38
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.	
Item dimensions (inches)	R820 small Standard base
Frame width	short legs: 21 long legs: 23
	short legs w/ casters: 22½ long legs w/ casters: 23½
Seat height above floor	short legs: 9½–12½ long legs: 18½–21½
	short legs w/ casters: 13½–16½ long legs w/ casters: 18½–21½
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7–12
Seat width with hip guides without hip guides	7–9 12
Seat depth with mini kit	8–12 7–11
Armrest height above seat	5–7½
Trunk support width	5½–11½
Backrest height	12½–15½
Headrest height above seat	14½–21
Max. working load (lbs)	75



Small hi/lo base
R830 Rifton Activity Chair

User dimensions (inches)	
Height with mini kit	32–48 28–38
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.	
Item dimensions (inches)	R830 small Hi/lo base
Frame width	26
Seat height above floor	10–25
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7–12
Seat width with hip guides without hip guides	7–9 12
Seat depth with mini kit	8–12 7–11
Armrest height above seat	5–7½
Trunk support width	5½–11½
Backrest height	12½–15½
Headrest height above seat	14½–21
Max. working load (lbs)	75

Room to grow

Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"




Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"




Medium standard base
R840 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
Frame width	short legs: 23¾
	long legs: 24½
Seat height above floor	short legs w/ casters: 24½
	long legs w/ casters: 25½
Seat height above footboard	short legs: 12½–16½
	long legs: 19–23
Seat angle (tilt-in-space)–degrees	short legs w/ casters: 16–20
	long legs w/ casters: 19–23
Backrest angle–degrees	15° forward, 15° back
Footboard knee angle–degrees	5° forward, 20° back
Seat height above footboard	45°–110°
Seat width with hip guides	7¾–15¼
without hip guides	8½–11½
Seat depth	14
Armrest height above seat	11–16
Distance between laterals	7–10½
Backrest height	6½–12
Headrest height above seat	15½–19
Max. working load (lbs)	17½–24½
	150



Medium hi/lo base
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides	8½–11½
without hip guides	14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lbs)	150

Room to grow

Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"




Large chair at its largest dimensions

- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"




Large standard base
R860 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	R860 large Standard base
Frame width	short legs: 26¾ long legs: 26¾
	short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24
	short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lbs)	250



Large hi/lo base
R870 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	R870 large Hi/lo base
Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lbs)	250