

Healthcare
furniture,
enhancing
lives.

Caremed Alrick.

8000 Series Companion Bed

Two patients can be cared for in comfort and dignity with the Alrick 8000 Companion Series bed. It features electric hi-low, electric knee brake, and a 300kg working load. The actuators are heavy duty, ensuring efficient, reliable function, and are warranted for five years. The steel bed frame is warranted for ten years. Available in both Double and Queen size.



Configuration Possibilities




OPTIONAL Upholstered Surround

Accessories Available




360HR Horizontal Rails **380VR** Vertical Rails **400HR** Deluxe Horizontal Rails **CLP400** Clamp on Support Rail



ASHP Self Help Pole & Triangle **AIVP** I.V. Pole **TBR2000** T-Bar Wall Protection Buffer


Standard Functions



Electric Hi-Low



Electric Backrest



Electric Kneebend

Standard features also include;

- Manual vascular support
- 125mm braking and directional castors
- Battery backup
- APO Grey powdercoat colour
- Fini Blonde laminate head and footboards

Specifications

Overall dimensions (mm)	2060 x 1450 Double
	2060 x 1580 Queen
Sleeping surface dimensions (mm)	1990 x 1370 Double
	1990 x 1476 Queen
SWL	300kg
NET Bed weight	150kg
Low Position	370mm
High Position	700mm
Underbed clearance	160mm
Angles (°)	Back 75
	Knee 55
Warranty	Frame 10 yr
	Electric 5yr
	Castors 1yr
	Handset 1 yr
Optional features	Underbed light
	Backlit Handset
	Central locking
	End board options
Product codes	Double 8000D
	Queen 8000Q

** LOAD CAPACITY HAS BEEN TESTED IN ACCORDANCE WITH AS/NZS 3200.2.38:2007 **