

Healthcare
furniture,
enhancing
lives.

Caremed Alrick.

Alrick EN7
Endless Bed Series

'True floorline bed
with no towers and
endless possibilities'



The EN7 series bed is a variant of the EN9000 manufactured overseas. It is the first ever universal bed that suits low care, high care and bariatric care, plus the innovative three-in-one bed. The EN7 uses Linak motors and has a low height of 120mm and has no towers!

Configuration Possibilities



Features

Electric Hi/Low	Auto regression backrest
Electric kneebend	75mm Twin wheel castors
Electric forward & reverse Trendelenburg	Standard head and foot boards (Fini Blonde laminate)

Specifications

Overall dimensions (mm)	2080 x 920 Single	
	2080 x 1090 King single	
Sleeping surface dimensions (mm)	2000 x 890 Single	
	2000 x 1070 King single	
SWL	200kg/300kg	
NET Bed weight	110kg	
Low Position	120mm	
High Position	725mm	
Underbed clearance	Clearview	
	Angles (°)	Back 62
		Knee 35
		Trendelenburg 15
	Reverse 15	
Warranty	Frame	10 yr
	Electrics	5 yr
	Handset & Castors	1 yr

Optional features

Underbed light
Battery backup
Dual braking
End board options

- 01 Electric forward & reverse Trendelenburg
- 02 Electric Kneebend
- 03 Auto regression backrest
- 04 Standard head & foot boards
- 05 10 year frame & 5 year electrical warranty
- 06 Electric Hi-Low
- 07 75mm Twin wheel castors

