

Healthcare  
furniture,  
enhancing  
lives.

# Caremed Alrick.

Alrick EN9  
Endless Bed

A bed with Endless  
possibilities




We are proud to introduce the first ever universal bed that suits low care, high care and bariatric care. Plus, this innovative three-in-one bed has a low height of 100mm and NO TOWERS.

From care of frail elderly or dementia patients, through to acute and bariatric settings, the patent pending Australian made EN9000 Endless Bed has you covered. One bed, endless possibilities!


## Features




Australian Made



Unique dual safe working load—Simply swap handsets




Low height: 100mm  
High height: 725mm



No towers. Up to 400mm shorter than other floor level beds



Universal solution for Low, High & Bariatric care




Chair position at the push of a button



10 year frame & 5 year linak warranty

## Configuration Possibilities



## Specifications

Overall dimensions (mm)	2080 x 920mm Single
	2080 x 1090mm King Single
	2080 x 1390mm Double
Sleeping surface dimensions (mm)	2000 x 890mm Single
	2000 x 1070mm King Single
	2000 x 1370mm Double
SWL	200kg/300kg
NET Bed weight	110kg
Low Position	100mm
High Position	725mm
Underbed clearance	Clearview
Angles (°)	Back 62
	Knee 35
	Trendelenberg 15
	Reverse 15
Warranty	Frame 10 yr
	Electrics 5 yr
	Handset & Castors 1 yr
Optional features	Underbed light
	Battery backup
	Dual braking
	End board options
Product codes	Single EN9000
	King Single EN9000KS
	Double EN9000D