


Healthcare
furniture,
enhancing
lives.

Caremed **Alrick.**

M9040IMP Mattress

Four Way Turn Mattress



This Alrick 4 Way Turn Mattress has a profiled surface to aid free flow of air and reduce surface tension. It has a two way stretch vapour permeable cover with sewn seams and zipped closure and is profiled on both sides. The M9040IMP is also available in a variety of sizes, Single & King Single.

Specifications

Overall dimensions	Single	1960 x 870 x 150mm
	King Single	1050 x 1960 x 150mm
Product codes	Single	M9040IMP
	King Single	M9040KSIMP
SWL	30-213kg	
Warranty	2 yr	

Features

- ✓ Profiled foam surface
- ✓ 2 Way stretch cover
- ✓ Vapour permeable cover
- ✓ Sewn Seams
- ✓ Zipped closure, profiled on both sides
- ✓ Risk category - Very High
- ✓ Flame resistant cover and foam
- ✓ 4 way turn

