

Healthcare
furniture,
enhancing
lives.

Caremed Alrick.


M9080IMP
Premium Pressure Care Mattress



The Alrick Premium Pressure Care Mattress is renowned for the fact that it reduces cleaning time which increases cover longevity. The Covers can be cleaned in-situ for easy care or laundered. It has been clinically tested, and the finest quality foams are used to provide excellent pressure distribution and improve patient comfort. It also aids movement on profiling bed frames and reduces shear and friction for ease of manual handling.

Specifications		
Overall Dimensions	Single	L 1200mm W 870mm H 150mm
	King Single	L 1050mm W 1960mm H 150mm
Weight Capacity	213kg	
Product Codes	Single	M9080IMP
	King Single	M9080KSIMP
Warranty	2 yr	

Features	
Profiled CMHR foam surface	High frequency welded cover
Machine-hinged profile design	Flame resistant cover & foam
Risk category - Very High	Vapour permeable, waterproof
Mattress cover material PU coated fabric, multi-stretch	Welded seams double sided welded flap



CMHR Foam



Welded Seams



Profiling Interior Hinges