PivotFit[™] MD Series Shoulder Harness

Promotes upright posture through shoulder control for improved breathing, head control and stability. Patented, body-contouring design distributes pressure across a 25% greater surface area than a traditional H-style harness, while the Swivel Buckle allows pads to pivot into place, easily accommodating asymmetry.



reduce choking risk.

Style Options



Dynamic Style Patented, dual-zone stretch enhances comfort by allowing moderate shoulder movement without reducing support at the chest. Ideal for helping to manage spastic tone.



Standard Style Fixed, non-stretch style reduces movement and has an easy-toclean CORDURA® nylon covering. Effective for severe kyphosis, spinal curvature or trunk rotation.

Part Numbers

Includes four end-fittings. Choose between high-strength Plastic Flat and Cinch-Mounts[™] and add suffix -B5 or -B6 to the end of the part number.

Size	PivotFit™ Dynamic	PivotFit™ Standard
XXS	SH295XXS	SH290XXS
XS	SH295XS	SH290XS
S	SH295S	SH290S
М	SH295M	SH290M
L	SH295L	SH290L
XL	SH295XL	SH290XL

End-Fittings

Plastic Flat-Mount[™]



Plastic Cinch-Mount™



Visit us at **bodypoint.com** 206.405.4555 or 800.547.5716 | sales@bodypoint.com | Fax: 206.405.4556 © 2020 Bodypoint. All Rights Reserved. BMM326 2020.3

