

TOTAL LENGTH

S38	127cm (50")
M40	140cm (55")
M46	152cm (60")

PADED LENGTH

S38	38cm (15")
M40	40cm (15-3/4")
M46	46cm (18")

CLOSED CELL FOAM PAD
S = 7mm THICK
M = 11mm THICK

WEBBING WIDTH

S38	2.5cm (1")
M40 & M46	4cm (1-1/2")

PAD WIDTH

S38	5cm (2")
M40 & M46	5.5cm (2-1/4")

ADJUSTMENT STRAP RANGE

S38	5cm (2")
M40	6.5cm (2-1/2")
M46	7.5cm (3")

SECONDARY STRAP LENGTH

S38	46cm (18")
M40	51cm (20")
M46	58.5cm (23")

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:




PUSH-BUTTON

SIZE	PART NUMBER
S38	N/A
M40	HB409-M40
M46	HB409-M46


Sub-ASIS COMPATIBLE BELT

M46	LB409-M46
-----	-----------



PLASTIC SIDE-RELEASE

SIZE	PART NUMBER
S38	HB419-S38
M40	HB419-M40
M46	HB419-M46



REHAB LATCH

SIZE	PART NUMBER
S38	N/A
M40	HB429-M40
M46	HB429-M46

Sub-ASIS COMPATIBLE BELT

M46	LB429-M46
-----	-----------



X-SMALL METAL PUSH-BUTTON

SIZE	PART NUMBER
S38	HB439-S38
M40	N/A
M46	N/A

-A1 + **FLAT-MOUNT™**
END FITTING PART NUMBERS
S (1") HKA1-R-1
M (1-1/2") HKA1-G-1
L (2") HKA1-B-1

-B1 +

-A2 + **CINCH-MOUNT™**
END FITTING PART NUMBERS
S (1") HKA2-R-1
M (1-1/2") HKA2-G-1
L (2") HKA2-B-1

-B2

note

description

HIP BELT, 4-POINT, REAR-PULL, TECHNICAL SPECS

drawing no.

D16-1116-32

date originated

11/16/2016

drawn by

D. DOUGAL

rev.

B