

Healthcare  
furniture,  
enhancing  
lives.

Caremed  
**Alrick.**

3004 Paediatric Changeable

Electric foot control operation



The Alrick 3004 Series Single Tower Paediatric Changetable is ideal for children and small adults up to 90kg. The foot control operation provides flexible and easy to control height adjustment, while the single tower lifting column conceals electrics within the column for safe and smooth hi-low functioning.

Specifications		
Height	Minimum	480mm
	Maximum	880mm
Overall dimensions	Width	700mm
	Length	1400mm
SWL	90kg	
Product code	3004	
Warranty	Frame	10 yr
	Electric	5 yr

Features



Foot Control

Electric Hi/Low	125mm Braking Castors
Blue Padded Vinyl Top	Foot Control Operation



Australian Made



10 year frame



5 year Linak warranty

Options



LAP/B Lap Belt



HB61000-00102 Hand Control



380VRKCT Side Rails