



# momo therapy bike.

## Invoice address:

## Dispatch address:

See invoice address

## Order

## Cost estimate

Customer No.: \_\_\_\_\_

Orderer: \_\_\_\_\_

Date: \_\_\_\_\_

Signature: \_\_\_\_\_

## Commission:

## The therapy bike for more balance.


- Aluminium frame with extra-low access point
- Rim brake with separate parking brake
- Angle adjustable handlebar stem
- Steering damper for stabilisation when moving straight on
- Drum brake in front wheel
- Stabiliser frame with dynamic supports and adjustable spring strength
- Incl. luggage carrier

**HMV-No.: 22.51.01.0014**






### Dimensions - basic model


		12"	16"	20"	24"	26"
<b>A*</b>	Pedal to upper edge of saddle	30 - 38 cm	42 - 57.5 cm	50 - 64.5 cm	59 - 76.5 cm	71.5 - 87.5 cm
<b>A**</b>	Pedal to upper edge of saddle	33 - 41 cm	45 - 60 cm	52 - 68 cm	61 - 84 cm	73 - 95 cm
<b>B</b>	Saddle support tube to front tube	33 cm	39 cm	45 cm	53 cm	58 cm
<b>C</b>	Wheelbase	64 cm	76 cm	89 cm	104 cm	112 cm
<b>D</b>	Total length	99 cm	122 cm	148 cm	168 cm	178 cm
<b>E</b>	Access point height	16 cm	18 cm	20 cm	22 cm	23 cm
<b>F</b>	Total width	43 cm	52 cm	62 cm	71 cm	76 cm
	Max. load	35 kg	50 kg	65 kg	100 kg	100 kg
	Weight	10.5 kg	13 kg	15 kg	18 kg	20 kg



\* Dimensions with standard saddle post / \*\* Dimensions with T-saddle post




Basic model									
momo therapy bike. 12"		momo therapy bike. 16"		momo therapy bike. 20"		momo therapy bike. 24"		momo therapy bike. 26"	
Basic features, without drum brake		Basic features		Basic features		Basic features, additional lighting system according to StVZO		Basic features, additional lighting system according to StVZO	
Price / piece		Price / piece		Price / piece		Price / piece		Price / piece	
<b>38 01 000</b> 1.018 €		<b>38 02 000</b> 1.142 €		<b>38 03 000</b> 1.245 €		<b>38 04 000</b> 1.490 €		<b>38 05 000</b> 1.490 €	

Frame colours									
Colour selection for 12"			Colour selection for 16"+20"			Colour selection for 24"+26"			
									
Price			Price			Price			
Telemagenta (Pink) <b>(92 00 012)</b>		incl.	Melon yellow <b>(92 00 006)</b>		incl.	Anthracit metallic <b>(92 00 014)</b>		incl.	
Sky blue <b>(92 00 013)</b>		incl.	Marine blue <b>(92 00 008)</b>		incl.	Pure white <b>(92 00 100)</b>		incl.	
Pure white <b>(92 00 100)</b>		incl.	Pure white <b>(92 00 100)</b>		incl.				

Configuration possibilities - drive			
	Art.-No.	Drive	Upcharge / piece
	<b>38 00 001</b>	Rigid sprocket (Standard configuration, if nothing else stated)	incl.
	<b>38 00 101</b>	Freewheel brake hub (with back pedal brake)	incl.
	<b>38 00 002</b>	Extra light transmission (only for above drive types)	incl.
	<b>38 00 003</b>	3-gear freewheel brake hub (with back pedal brake)	199 €
	<b>38 00 007</b>	3-gear freewheel brake hub with reverse option	286 €
	<b>38 00 005</b>	7-Gang-Freewheel brake hub (with back pedal brake)	330 €
	<b>38 00 008</b>	7-gear-freewheel hub with roll brake (Brake lever mounted on the handlebar, without back pedal brake)	475 €

Saddle											
<b>Standard saddle</b> (standard configuration, if nothing else stated)			<b>Gel saddle</b>			<b>Saddle seat with leg hollows</b>					
											
	Width	Length		Rear width	Front width	Length		Rear width	Front width	Rear length	Total length
Size 1	15cm	21cm									
Size 2	15cm	24cm	Size 1	19cm	4cm	24cm	Size 1	20.5cm	4cm	9.5cm	14.5cm
Size 3	18cm	26cm	Size 2	24cm	7cm	27cm	Size 2	24cm	4cm	9.5cm	14.5cm
Price / piece			Upcharge / piece			Upcharge / piece					
Size 1 <b>(37 01 001)</b>		incl.	Size 1 <b>(37 01 024)</b>		50 €	Size 1 <b>(37 01 003)</b>		191 €			
Size 2 <b>(37 02 001)</b>		incl.	Size 2 <b>(37 02 024)</b>		50 €	Size 2 <b>(37 02 003)</b>		191 €			
Size 3 <b>(37 03 001)</b>		incl.									

Saddle					Saddle supports										
<b>Unicycle saddle</b> – Banana shape (caution: min. inside leg length increases by 2 cm)					<b>Moped saddle</b> (caution: min. inside leg length increases by 2 cm)				<b>Standard saddle support</b> (standard configuration, if nothing else stated)			<b>T-saddle support</b> – for horizontal saddle adjustment			
															
Rear width	Middle width	Front width	Length		Rear width	Front width	Rear length	Total length		for	for		for		
11cm	6cm	8cm	25cm		26cm	9.5cm	12cm	25cm		Size 1	Size 1		Size 1		
										Size 2	Size 2		Size 2		
										Size 3	Size 3		Size 3		
Upcharge / piece					Upcharge / piece				Price / piece			Upcharge / piece			
<b>37 00 010</b> 93 €					<b>37 00 023</b> 139 €				Size 1 <b>(37 01 004)</b> incl.			Size 1 <b>(37 01 005)</b> 73 €			
									Size 2 <b>(37 02 004)</b> incl.			Size 2 <b>(37 02 005)</b> 73 €			
									Size 3 <b>(37 03 004)</b> incl.			Size 3 <b>(37 03 005)</b> 73 €			

Holding brackets				Headrest			
<b>Holding brackets incl. holder</b> – only for use together with the T-saddle support				<b>Headrest</b> – only for use together with a holding bracket			
	<b>max. Height*</b>	<b>Depth</b>	 		<b>Width</b>	<b>Height</b>	
<b>Size 1</b>	20 cm	12 cm					
<b>Size 2</b>	30 cm	12 cm					
<b>Size 3</b>	37 cm	12 cm					
<b>Size 4</b>	53 cm	12 cm					
<b>Tailored</b>	___ cm (Max. 85 cm)	___ cm (Max. 30 cm)			<b>Size 1</b>	20 cm	15 cm
* Max. height: measured from the standard saddle to the top edge of the respective support pad				<b>Size 2</b>	23 cm	18 cm	
				Price / set		Price / piece	
Size 1	<b>(37 01 007)</b>		114 €	Size 1	<b>(37 01 029)</b>		303 €
Size 2	<b>(37 02 007)</b>		114 €	Size 2	<b>(37 02 029)</b>		303 €
Size 3	<b>(37 03 007)</b>		114 €				
Size 4	<b>(37 04 007)</b>		114 €				
Tailored	<b>(37 09 007)</b>		240 €				

Torso supports				Pelvic supports				
<b>Dynamic back support pad</b> padded with strap and bracket (only for use together with a holding bracket)		<b>Width adjustable back support pad</b> padded. incl. bracket (only for use together with a holding bracket)		<b>Dynamic pelvic guide support pad</b> padded with strap and bracket (only for use together with a holding bracket)		<b>Width adjustable pelvic guide support pad</b> – padded. incl. bracket (only for use together with a holding bracket)		
								
	<b>Width</b>		<b>Width</b>		<b>Width</b>		<b>Width</b>	
<b>Size 1</b>	20 - 30 cm	<b>Size 1</b>	20 - 28 cm	<b>Size 1</b>	20 - 30 cm	<b>Size 1</b>	20 - 28 cm	
<b>Size 2</b>	25 - 35 cm	<b>Size 2</b>	25 - 33 cm	<b>Size 2</b>	25 - 35 cm	<b>Size 2</b>	25 - 33 cm	
<b>Size 3</b>	30 - 40 cm	<b>Size 3</b>	30 - 36 cm	<b>Size 3</b>	30 - 40 cm	<b>Size 3</b>	30 - 36 cm	
Price / piece		Price / piece		Price / piece		Price / piece		
Size 1	<b>(37 01 008)</b>	213 €	Size 1	<b>(37 01 025)</b>	286 €	Size 1	<b>(37 01 009)</b>	213 €
Size 2	<b>(37 02 008)</b>	213 €	Size 2	<b>(37 02 025)</b>	286 €	Size 2	<b>(37 02 009)</b>	213 €
Size 3	<b>(37 03 008)</b>	213 €	Size 3	<b>(37 03 025)</b>	286 €	Size 3	<b>(37 03 009)</b>	213 €
						Size 1	<b>(37 01 026)</b>	286 €
						Size 2	<b>(37 02 026)</b>	286 €
						Size 3	<b>(37 03 026)</b>	286 €

Harness system								
<b>Chest strap</b> – for width adjustable back support pad				<b>Positioning jacket</b> – for width adjustable back support pad, incl. Universal bracket				
								
	<b>Width</b>	<b>Length</b>		<b>Inside width</b>	<b>Total width</b>	<b>Side height</b>	<b>Total height</b>	
<b>Size 1</b>	5 cm	18.5 cm		<b>Size 1</b>	8 cm	25 cm	10 cm	28 cm
<b>Size 2</b>	6.5 cm	25.5 cm		<b>Size 2</b>	9 cm	30 cm	12 cm	30 cm
<b>Size 3</b>	7 cm	30.5 cm		<b>Size 3</b>	12 cm	34 cm	16 cm	35 cm
Price / piece				Price / piece				
Size 1	<b>(37 01 050)</b>		86 €	Size 1	<b>(37 01 051)</b>		213 €	
Size 2	<b>(37 02 050)</b>		86 €	Size 2	<b>(37 02 051)</b>		213 €	
Size 3	<b>(37 03 050)</b>		86 €	Size 3	<b>(37 03 051)</b>		213 €	

**Harness system**

T-shape groin harness – for width adjustable pelvic support pad					4-point pelvic harness – for width adjustable pelvic support pad						
	Width of seat	Front width	Length of seat	Total length			Width of outer parts	Width of centre section	Outer part length	Centre section length	
<b>Size 1</b>	24 cm	23 cm	15 cm	32 cm		<b>Size 1</b>	7 cm	7 cm	12 cm	12 cm	
<b>Size 2</b>	26 cm	26 cm	17 cm	34 cm		<b>Size 2</b>	9 cm	8,5 cm	13 cm	14 cm	
<b>Size 3</b>	28 cm	28 cm	19 cm	38 cm		<b>Size 3</b>	11,5 cm	11 cm	15 cm	16 cm	
Price / piece					Price / piece						
Size 1 ( <b>37 01 052</b> )	146 €				Size 1 ( <b>37 01 053</b> )	150 €					
Size 2 ( <b>37 02 052</b> )	146 €				Size 2 ( <b>37 02 053</b> )	150 €					
Size 3 ( <b>37 03 052</b> )	146 €				Size 3 ( <b>37 03 053</b> )	150 €					

**Accessories**

Universal bracket for harness used without support pads (width = 27,5 cm)		Hand positioner – for all handlebars (if no quantity was specified: standard = 1 piece)		Steering stop – adjustable up to direction lock (incl. for steering for attendant)		Handlebar extension						
	Front width		Menge		for		Length					
	<b>Size 0</b>							8,5 cm	<b>Size 1</b>	12"	<b>Size 1</b>	10 cm
	<b>Size 1</b>							10 cm	<b>Size 2</b>	16"–26"	<b>Size 2</b>	10 cm
	<b>Size 2</b>							12 cm				
<b>Size 3</b>	14 cm											
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece	Price / piece	Price / piece	Price / piece					
<b>37 02 055</b>	97 €	Size 0 ( <b>37 00 054</b> )	62 €	Size 1 ( <b>37 01 006</b> )	114 €	Size 1 ( <b>37 01 022</b> )	25 €					
		Size 1 ( <b>37 01 054</b> )	62 €	Size 2 ( <b>37 03 006</b> )	114 €	Size 2 ( <b>37 02 022</b> )	25 €					
		Size 2 ( <b>37 02 054</b> )	62 €									
		Size 3 ( <b>37 03 054</b> )	62 €									

**Handlebar**

Classic handlebar – black (standard configuration, if nothing else stated)		Touring handlebar		Round handlebar		Multifunctional handlebar suitable for 20" – 26"								
	Standard for	Width		Width	Depth									
	<b>Size 1</b>	12"						47 cm	<b>Size 1</b>	55 cm	17 cm	<b>Size 1</b>	61,5 cm	16,5 cm
	<b>Size 2</b>	16"						49 cm	<b>Size 2</b>	58 cm	17 cm	<b>Size 2</b>		
	<b>Size 3</b>	20"						52 cm						
<b>Size 4</b>	24"+26"	58 cm												
Price / piece		Upcharge / piece		Upcharge / piece		Upcharge / piece								
Size 1 ( <b>37 01 011</b> )	incl.	Size 1 ( <b>37 01 012</b> )	98 €	Size 1 ( <b>37 01 013</b> )	98 €	<b>37 02 014</b>	150 €							
Size 2 ( <b>37 02 011</b> )	incl.	Size 2 ( <b>37 02 012</b> )	98 €	Size 2 ( <b>37 02 013</b> )	98 €									
Size 3 ( <b>37 03 011</b> )	incl.			Size 3 ( <b>37 03 013</b> )	98 €									
Size 4 ( <b>37 04 011</b> )	incl.													

**Handlebar brake**

Round handlebar with brake lever ring			Round handlebar with brake function – the brake function is simply activated by pressing down the round handlebars				
	Width	Depth		Width	Depth		
	<b>Size 1</b>	40 cm		17 cm	<b>Size 1</b>	40 cm	17 cm
	<b>Size 2</b>	43 cm		21 cm	<b>Size 2</b>	43 cm	21 cm
	<b>Size 3</b>	48 cm		25 cm	<b>Size 3</b>	48 cm	25 cm
Price / piece			Price / piece				
Size 1 ( <b>37 01 028</b> )	276 €		Size 1 ( <b>37 01 031</b> )	276 €			
Size 2 ( <b>37 02 028</b> )	276 €		Size 2 ( <b>37 02 031</b> )	276 €			
Size 3 ( <b>37 03 028</b> )	276 €		Size 3 ( <b>37 03 031</b> )	276 €			

## Accessories

<b>Push bar</b> height adjustable and removable, with integrated twist lock 	<b>Parking brake for attendant</b> (can only be used in connection with a push bar) 	<b>Parking brake for attendant</b> incl. motor switch-off when using the e-drive (can only be used in conjunction with Push bar)	<b>Steering for attendant</b> via deflection bar, incl. steering stop (can only be used in connection with a push bar)
Price / piece	Price / piece	Price / piece	Price / piece
for 12" <b>(3801009)</b> 158 € for 16" <b>(3802009)</b> 158 € for 20"-26" <b>(3803009)</b> 158 €	<b>3800011</b> 105 €	<b>3801011</b> 105 €	for 12" <b>(3801010)</b> 398 € for 16" <b>(3802010)</b> 398 € for 20" <b>(3803010)</b> 398 €

## Crank

<b>Progressively adjustable crank reducer</b> without widening the pedal distance 	<b>Crank reducer</b> – minimal pedal distance widening, reduction in size by 2.5 or 5 cm respectively 	<b>Special crank for knee contracture</b> for 20" – 26" 
Price / pair	Price / pair	Price / piece
for 12"-20" <b>(3700016)</b> 156 € for 24"+26" <b>(3704016)</b> 156 €	for 16"+20" <b>(3702015)</b> 80 € for 24"+26" <b>(3705015)</b> 80 €	mounted on right <b>(3703015)</b> 332 € mounted on left <b>(3704015)</b> 332 €

## Pedals

## Foot cups


<b>Foot fixing pedals</b> – with racing hooks and compensation weights 	<b>Home exerciser pedals</b> 	<b>Foot cups</b> – width adjustable, with fixing straps 						
Depth	Width		Min. front width	Max. front width	Min. rear width	Max. rear width	Length	
13cm	12cm		Size 1	8.7 cm	10.4 cm	5.7 cm	7.4 cm	17.4 cm
			Size 2	9.5 cm	12 cm	6.7 cm	9.2 cm	20.1 cm
			Size 3	11.5 cm	14 cm	8.1 cm	10.6 cm	23.8 cm
Upcharge / pair	Upcharge / pair	Price / pair						
<b>3700013</b> 84 €	<b>3700014</b> 58 €	Size 1 <b>(3701017)</b> 197 € Size 2 <b>(3702017)</b> 197 € Size 3 <b>(3704017)</b> 197 €						

## Foot cups with leg guidance

Foot cups with leg guidance – width adjustable and with fixing straps								
	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height	
<b>Size 1</b>	8.7 cm	10.4 cm	5.7 cm	7.4 cm	17.4 cm	15 cm	18.5 cm	
<b>Size 2</b>	9.5 cm	12 cm	6.7 cm	9.2 cm	20.1 cm	18 cm	22 cm	
<b>Size 3</b>	11.5 cm	14 cm	8.1 cm	10.6 cm	23.8 cm	21 cm	26 cm	
<b>Size 4</b>	11.5 cm	14 cm	8.1 cm	10.6 cm	23.8 cm	28 cm	36 cm	
	Price / pair							
Size 1 <b>(3701018)</b>								334 €
Size 2 <b>(3702018)</b>								334 €
Size 3 <b>(3703018)</b>								334 €
Size 4 <b>(3704018)</b>								334 €

**Foot cups with dynamic leg guidance**

**Foot cups with dynamic leg guidance**  
width adjustable and with fixing straps




	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height
<b>Size 1</b>	8.7 cm	10.4 cm	5.7 cm	7.4 cm	17.4 cm	15 cm	18.5 cm
<b>Size 2</b>	9.5 cm	12 cm	6.7 cm	9.2 cm	20.1 cm	18 cm	22 cm
<b>Size 3</b>	12 cm	15 cm	8.8 cm	11.8 cm	24 cm	20.5 cm	25.5 cm

Price / pair

Size 1 <b>(37 01 035)</b>	437 €
Size 2 <b>(37 02 035)</b>	437 €
Size 3 <b>(37 03 035)</b>	437 €

**Configuration possibilities - lighting system**

	Art.-No.	Lighting system	for	Upcharge / set
	<b>38 01 021</b>	Upgrade with a lighting system with side runners and reflectors according to StVZO	16"+20"	170 €
	<b>38 02 021</b>	Upgrade with a lighting system with hub dynamo and reflectors according to StVZO (not possible in conjunction with an electric drive)	16"+20"	286 €
	<b>38 05 021</b>	Lighting system in conjunction with an electric drive (this lighting system is already included on 24" + 26"-tricycles, which are ordered with an electric drive)	16"+20"	170 €
	<b>38 03 021</b>	Lighting system with hub dynamo (not possible in conjunction with an electric drive)	24"+26"	122 €

**Accessories**
**Basket**


Price / piece

for 12"+16"	<b>(38 01 019)</b>	50 €
for 20"+24"	<b>(38 03 019)</b>	50 €
for 26"	<b>(38 05 019)</b>	50 €

**Quick connection of stabiliser frame** – for simple disassembly for transport

Price / set

for 12"	<b>(38 01 026)</b>	118 €
for 16"	<b>(38 02 026)</b>	118 €
for 20"	<b>(38 03 026)</b>	118 €
for 24"	<b>(38 04 026)</b>	118 €
for 26"	<b>(38 05 026)</b>	118 €


**Upgrade set**

**Stabiliser frame for third-party bicycles** – with dynamic supports and adjustable spring strength

Price / set

for 12"	<b>(38 01 101)</b>	606 €
for 16"	<b>(38 02 101)</b>	606 €
for 20"	<b>(38 03 101)</b>	606 €
for 24"	<b>(38 04 101)</b>	663 €
for 26"	<b>(38 05 101)</b>	663 €
Assembly at Schuchmann <b>(38 00 025)</b>		Net 145 €

**E - drive**
**DirectPower E-drive from Heinzmann**

The new **momo therapy bike** is equipped with the DirectPower E-drive from Heinzmann. This allows us to offer you a wear-free and silent direct drive. The motor is designed to be brushless and without gears. This minimises maintenance and wear and eliminates motor noise. The recuperative system automatically recharges when driving downhill or braking. This means that the range can be increased by up to 15 %. The existing force of the user is measured by a force sensor integrated in the bottom bracket, and the assistance from the electric drive is adjusted.

Choose the suitable E-drive for the **momo therapy bike** from the following configurations. You can choose whether the drive should be equipped **with or without start-up assistance**. In addition, you can choose between a **beginner version** and an **advanced version**. The difference here is the **start-up** and the **top speed** of the various modes **Eco**, **Standard** and **Power**, which can be selected via the display. An analysis cable can be used to change the selected version at a later date.



The e-drive excludes a rigid gear rim and brings with it an additional weight of approx. 10 kg.

**E - drive configuration**
**DirectPower E-drive from Heinzmann – without start-up assistance**

for size	Version		Start-up handling / top speed			Price / piece
			Eco	Standard	Power	
16"	Beginner	<b>(38 02 085)</b>	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	2.540 €
	Advanced	<b>(38 02 086)</b>	medium / 6 km/h	strong / 8 km/h	strong / 10 km/h	2.540 €
20"	Beginner	<b>(38 03 085)</b>	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	2.540 €
	Advanced	<b>(38 03 086)</b>	medium / 6 km/h	strong / 8 km/h	strong / 12 km/h	2.540 €
24"	Beginner	<b>(38 04 085)</b>	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	2.540 €
	Advanced	<b>(38 04 086)</b>	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	2.540 €
26"	Beginner	<b>(38 05 085)</b>	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	2.540 €
	Advanced	<b>(38 05 086)</b>	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	2.540 €

**E - drive configuration**
**DirectPower E-drive from Heinzmann – with start-up assistance**

Select whether the start-up assistance should be activated by pressing the button (**Fig. 1**) or by twist handle (**Fig. 2**).



for size	Version	Start-up assistance	Start-up handling / top speed			Price / piece	
			Eco	Standard	Power		
16"	Beginner with twist handle	<b>(38 02 087)</b>	4 km/h	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	2.643 €
	Beginner with button	<b>(38 02 088)</b>		gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	2.643 €
	Advanced with twist handle	<b>(38 02 089)</b>	6 km/h	medium / 6 km/h	strong / 8 km/h	strong / 10 km/h	2.643 €
	Advanced with button	<b>(38 02 090)</b>		medium / 6 km/h	strong / 8 km/h	strong / 10 km/h	2.643 €
20"	Beginner with twist handle	<b>(38 03 087)</b>	4 km/h	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	2.643 €
	Beginner with button	<b>(38 03 088)</b>		gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	2.643 €
	Advanced with twist handle	<b>(38 03 089)</b>	6 km/h	medium / 6 km/h	strong / 8 km/h	strong / 12 km/h	2.643 €
	Advanced with button	<b>(38 03 090)</b>		medium / 6 km/h	strong / 8 km/h	strong / 12 km/h	2.643 €
24"	Beginner with twist handle	<b>(38 04 087)</b>	6 km/h	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	2.643 €
	Beginner with button	<b>(38 04 088)</b>		gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	2.643 €
	Advanced with twist handle	<b>(38 04 089)</b>	6 km/h	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	2.643 €
	Advanced with button	<b>(38 04 090)</b>		gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	2.643 €
26"	Beginner with twist handle	<b>(38 05 087)</b>	6 km/h	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	2.643 €
	Beginner with button	<b>(38 05 088)</b>		gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	2.643 €
	Advanced with twist handle	<b>(38 05 089)</b>	6 km/h	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	2.643 €
	Advanced with button	<b>(38 05 090)</b>		gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	2.643 €

**Software**
**Subsequent conversion/adjustment of the software on the delivered product by our sales force**

Net price (not discountable) / piece

**38 00 091**

134 €