

# momo tricycle.

## Let's go for a little ride.

GMFCS  
Levels I-V

kg  
max. 120




30-98.5 cm

16-23 cm

**I am your cool tricycle.** I can be there for you whenever you want to take a ride around your streets or a lake. It's fun, but we also train and mobilise your joints, musculature and your circulation. You can also select me in very different versions, saddle sizes, handlebar types and colours. Your age doesn't matter, because I am available in lots of different sizes: for children, adolescents and adults.





# momo tricycle.



**Support**  
for keeping the upper body upright.

**Saddle**  
available in different shapes and sizes.

**Handlebars**  
available in many different designs.

**Optimised geometry**  
for easy transfer, a good sitting position and high manoeuvrability.

**High-quality chain guide**  
runs extremely quietly.

**Individual drive**  
from a rigid drive right up to 7-gear freewheel brake hub.

**Low access level as standard feature**  
permits the rider to climb on conveniently and also allows independent transfer.

**Frame**  
made from lightweight, highly robust aluminium.

You can find a diverse range of accessories for optimum support & individual equipment from page 122.

**You have the choice! Between lots of handlebars, saddles, foot pans, head supports and different drives!**



**Frame colours 12"**

- Telemagenta
- Sky blue
- Pure white

**Frame colours 16" and 20"**

- Melon yellow
- Navy blue
- Pure white

**Frame colours 24", 26" and 26" (XL)**

- Anthracite metallic
- Pure white



Data	12"	16"	20"	24"	26"	26" (XL)
Pedal to top edge of saddle*	30 - 38 cm	42 - 57.5 cm	50 - 64.5 cm	56.5 - 73.5 cm	66.5 - 83 cm	74 - 90.5 cm
Pedal to top edge of saddle**	33 - 41 cm	45 - 60 cm	52 - 68 cm	58 - 81.5 cm	67.5 - 91 cm	75 - 98.5 cm
Saddle support tube to front tube	33 cm	39 cm	45 cm	52 cm	58 cm	58 cm
Wheelbase	68 cm	80 cm	91 cm	104 cm	112 cm	112 cm
Access point height	16 cm	18 cm	20 cm	22 cm	23 cm	23 cm
Full length	102 cm	120 cm	140 cm	165 cm	178 cm	178 cm
Full width	58 cm	65 cm	69 cm	75 cm	75 cm	75 cm
max. load	40 kg	60 kg	80 kg	100 kg	120 kg	120 kg
Tare weight	12 kg	14.5 kg	17 kg	19.5 kg	21.5 kg	21.7 kg
HMV No.	22.51.02.0053	22.51.02.0054	22.51.02.0055	22.51.02.0051	22.51.02.0052	22.51.02.0052

\* Dimensions with standard saddle post / \*\* Dimensions with T-saddle post



You can find further information on accessories & harnesses on pages 132 – 135 and in the price list.

