


Technical Specifications Floorstand



Model





FS
Floor Stand



FS
Floor Stand
Quickpack

Technical Specifications



Size & Weight

Floorstand #17.1050

min. Height: 103 cm / 40"
max. Height: 173 cm / 68"
Width: 77 cm / 30"
Length: 79 cm / 31"
Weight: 7,5 kg / 16.5 lbs.

Floorstand Quickpack #17.1052

min. Height: 108 cm / 42"
max. Height: 178 cm / 70"
Width: 77 cm / 30"
Length: 79 cm / 31"
Weight: 8,5 kg / 18.7 lbs.

Max. Load capability
6 kg / 13 lbs.

*optional HD-kit extends weight
limit to 10 kg / 22 lbs.*

Made by REHAdapt Engineering GmbH & Co. KG

Heinrich-Hertz-Str. 104 • 34123 Kassel • Germany

Tel.: +49 561 51 09 04 22 • Fax: +49 561 51 09 04 21 • office@REHAdapt.com • www.REHAdapt.com